ENGLISH

Indication:

The product is indicated for conditioning the internal surface of acid-sensitive dental ceramics, in order to promote adequate adhesion to resin systems. It is also recommended for the external repair of acid-sensitive ceramic restorations using composite resinous materials.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who require acid-sensitive ceramic restorations or the repair of previously installed ceramic restorations, regardless of age, gender, or systemic condition, as long as there is a clinical indication for restorative procedures with resin adhesive systems.

Fundamentals of its operation and action:

The hydrofluoric acid gel conditioner is a 10% gel used for conditioning the internal and external surfaces of dental ceramics (conditionable ceramics). Hydrofluoric acid conditioning acts on the glass phase of ceramics that are acid-sensitive, promoting a penetration pattern on the surface, making it micro-retentive and generating microporosities. This leads to greater retention of ceramics, creating better adhesion with resinous materials, such as adhesives and cements.

Contraindications:

- Do not use to condition dental tissues (enamel, dentin, cementum and bone).
- Not for use on ceramics that are **not** acid-sensitive.
- Contraindicated in patients with a history of hypersensitivity to the components of the formulation.
- Do not use on a restorative surface that does not require acid etching.

Instructions for use:

The external surface conditioning (for composite repair) or internal (for cementation), depends on the type of ceramic, so it is recommended to consult the manufacturer's indication of the ceramic, we suggest the following times for internal conditioning:

- Feldspathic Porcelain: condition for 120 seconds;
- Reinforced vitreous ceramic (leucite): condition for 60 seconds;
- Vitreous ceramic (lithium disilicate): condition for 20 seconds;

And for external conditioning (repair):

- 60 to 90 seconds.

After homogeneous application of the product for surface conditioning (internal or external), proceed as follows for composite repair or cementation of prosthetic parts:

1. Initially remove all the acid with a suction tip;
2. Wash thoroughly and dry the conditioned surface;
3. Proceed with silanization according to the manufacturer's instructions;
4. Apply the adhesive system, resin cement or composite resin according to the manufacturer's instructions;

- Only in the case of repairs with composite resin may it be necessary to make a bevel of about 2mm in the fractured area with a diamond tip prior to the application of the acid.

ÁCIDO FLUORÍDRICO 10%
HYDROFLUORIC ACID 10%
ÁCIDO FLUORHÍDRICO 10%

REV.01

Form of presentation

- Product presented in the form of a gel;
- Syringe with 2.5g.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As a chemical, hydrofluoric acid should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

- Inspect if the product is in perfect condition for use.
- Remaining on the prosthetic specimen for more than specified times may cause excessive conditioning.
- The dentist and the patient must wear protective goggles. The dentist should also wear protective gloves.
- When the use is intra-oral (in cases of repair with composite resin) it is necessary to use insulation with rubber dam. If it is not possible to use it, make sure that the relative insulation is secure to apply the product.
- Because it contains a strongly acidic product, it should not come into contact with eyes, skin and oral tissues. In case of accidental contact, rinse the affected area with plenty of water.
- If irritations persist, seek medical advice.
- Do not use the product after the expiration date.
- Keep out of reach of children.
- Do not reuse the syringe.
- Non-injectable.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Useful life:

Single-use product, with a shelf life restricted to the moment of clinical application. After being removed from the original packaging and used, any remaining quantity should be discarded immediately.

Conservation and storage:

- Always keep the package tightly closed.
- Store the product in a cool, dry place.
- Do not exceed a temperature of 37 °C.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

Additional information:

This review cancels and replaces the previous revision
3R Indústria e Comércio Ltda.
Rua Ptolomeu, 290 – Vila Socorro
São Paulo - SP – Brasil - CEP: 04762-040

MS: 80676920037

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ÁCIDO FLUORÍDRICO 10%
HYDROFLUORIC ACID 10%
ÁCIDO FLUORHÍDRICO 10%

REV.01

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Hydrofluoric Acid

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Maker	Indicates the manufacturer of the health product as defined in Directive 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.

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