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ENGLISH

Indication:

It is indicated for the conditioning of enamel and/or dentin prior to the application of adhesives and consequently of restorative materials to the tooth such as composite resin, resin cement and resin sealants.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who will undergo restorative procedures that require prior acid etching, including direct resin composite restorations, adhesive cementations (such as fiberglass pins, crowns, veneers, inlays/onlays), and application of resin sealants. It covers adolescents and adults, as well as children when there is a specific clinical indication for adhesive procedures.

Fundamentals of its operation and action:

It is a thixotropic gel based on water and 37% phosphoric acid and low viscosity. It has a blue color for easy viewing and control when applied. It is fully soluble in water for easy removal after application. It contains phosphoric acid content suitable for a perfect chemical attack on the surface, promoting demineralization and creating microporosities that will later be filled by the adhesive, in addition to the removal of the "smear layer", this step helps in the adhesion of restorative materials to the tooth.

Contraindications:

- Contraindicated in patients with a history of hypersensitivity to the components of the formulation.
- Do not use on a surface that will be restored that does not require acid etching.

Instructions for Use:

- Perform isolation of the operative field;
- Perform prophylaxis with water and pumice stone;
- Apply the acid conditioner on enamel and/or dentin, remembering that in dentin you should wait 15 seconds and enamel 30 seconds;
- Apply a jet of abundant water until the product is completely removed;
- Dry the enamel with an air blast and the dentin with absorbent paper or cotton;
- Poor hygiene of the area, lack of absolute isolation, variable degree of fluorosis and imperfect amelogenesis can cause failure of partial or total conditioning.

Form of presentation

- Product presented in the form of a gel;
- Syringe with 2.5g.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As a chemical, phosphoric acid should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

- Because it contains a strongly acidic product, it should not come into contact with eyes, skin and oral tissues,
- Protect neighboring teeth with appropriate materials;
- In case of accidental contact, wash the affected area with plenty of water; if irritations persist, seek medical advice;
- Do not reuse the syringe;
- Do not use the product after the expiration date;
- Keep out of reach of children;
- Product for the exclusive use of the professional;
- During operation, the professional must use all personal protective equipment, including protective goggles;

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REV.01

- Wear protective glasses on the patient;
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer;

Disposal:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

Single-use product, with a service life restricted to the moment of clinical application. After being removed from the original packaging and used, any remaining quantity should be discarded immediately.

Conservation and storage:

- Keep the packaging tightly closed at all times
- Do not exceed a temperature of 37 °C
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Phosphoric acid

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).

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REV.01

	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).
---	--	--

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