AQUALIZE

REV.01

ENGLISH

Indication:

Aqualize is indicated for use in conjunction with ultrasound equipment in procedures that require cooling in the insert. With a 1-liter reservoir, the Aqualize ensures continuous cooling throughout the procedure, eliminating the need to use the office's water flow. Thus, the equipment guarantees not only efficiency and safety, but also versatility.

Fundamentals of its operation and action:

Aqualize is a portable automatic water supply system, which allows the use of ultrasound equipment without relying on water from the dental office. Proper refrigeration while using ultrasound is critical to prevent overheating of dental tissues and associated damage.

EQUIPMENT INSTALLATION

Open the package and follow the instructions in the following illustration:

- Connect the device to the Ultrasound (Figure 1).
- Pour water into the bottle and tighten the bottle cap securely.
- Plug the power supply into the machine.
- Connect the power supply to the power outlet.
- Turn on the power button to start the device
- Turn off the equipment after finishing the operation.

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Instructions for use:

- Keep the product in a place protected from humidity and excessive heat, including the incision of sun rays on it.
- During product accommodation and operation, check that the environmental conditions are within the specified ranges as described below:
 - Ambient temperature: +10°C to +40°C.
 - Relative humidity: 0% to 80%.
 - Atmospheric pressure: 70 kPa – 106 kPa.

Cleaning Requirements:

- Use gauze or tissues appropriate for cleaning in the office.
- Use a non-corrosive disinfectant with a neutral PH value to clean the bottle and the pipe of water or 70% alcohol.

WARNING AND PRECAUTIONS TO BE TAKEN

- Ensure that all pipes are connected accurately and firmly before operating.
- Turn off the power source when you put in or change the water.
- Do not exceed the maximum volume of water in the bottle when pouring water.
- When the model is working, do not pull the tubes hard to avoid damaging the air tube and water tube.

AQUALIZE

REV.01

- Since there is air pressure in the bottle, take out the air tubing and water tubing before opening the bottle cap and injecting water.
- Please tighten the bottle cap tightly after changing the water and turn on the power switch.
- Do not remove the water pipe after completing pressurization to avoid splashing water from the bottle and causing secondary hazards.
- Please use our company's special adapter.
- Repair of the equipment must be performed by the manufacturer or authorized dealer.
- Do not use in the presence of flammable anesthetic gases or other flammable materials, such as some types of cleaning products. Do not start the equipment in areas where there may be a danger of explosion.
- Equipment not suitable for use in oxygen-rich environments.

ATTENTION: The equipment has an external power supply and should be treated as part of the equipment, do not use any source other than the one provided with the equipment, if you have any questions, contact your Microdont representative.

Conservation and storage:

- Keep the product in its original packaging and in a place protected from humidity and excessive heat.
- During transport, avoid vibrations, drops and impacts directly on the packaging.
- During transport and storage, ensure that environmental conditions are within the specified ranges as described below:
 - Ambient temperature for transport and storage: -20 °C to +55 °C.
 - Relative humidity for transport and storage 0% to 80%.
 - Atmospheric pressure: 50 kPa to 106 kPa

TECHNICAL SPECIFICATIONS

Electric

- Power Supply Voltage Input:
 - 100-240Vac/50/60Hz.
- Control Unit Voltage Input:
 - 5.0Vc.c/1.0A.

General

- Main unit weight: ~272g.
- Water volume capacity: 1L.
- Main Unit (Applicator) Dimensions: ~78x280mm (Cylindrical).
- IP protection level: Power supply - IPX0, Handpiece IPX0, Charging dock - IP00.
- Power supply weight: ~0.5Kg.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

AQUALIZE

REV.01

Simbologia:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Care	Indicates the need for the user to consult the instructions for use for important preventive information, such as warnings and precautions that cannot, for various reasons, be presented in the product itself for health.
	See instructions for use	Indicates the need for the user to refer to the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured.
	Keep away from sunlight	Indicates a health product that needs protection against light sources.
	Manufacturer	Indicates the Manufacturer of the product.
	Batch Code	Indicates the manufacturer's batch code so that the batch or matchup can be identified.
	Catalog number	Indicates the manufacturer's reference number so that the health product can be identified.
	Serial number	Indicates the serial number used to identify a specific device.

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