BICARBONATO DE SÓDIO
SODIUM BICARBONATE
BICARBONATO DE SODIO

REV.01

ENGLISH

Indication:

Sodium bicarbonate is indicated for dental prophylaxis procedures, assisting in the mechanical removal of bacterial plaque and extrinsic stains through application with air-polishing (air-abrasion) equipment.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who require removal of biofilm and extrinsic stains through mechanical prophylaxis, including individuals with good periodontal condition and no contraindications to the use of sodium bicarbonate air polishing.

Fundamentals of its operation and action:

Dental sodium bicarbonate acts as an abrasive agent in air-polishing equipment, promoting the controlled removal of plaque and extrinsic stains. Its mechanical action contributes to the cleaning of tooth surfaces and, in some cases, to the polishing of enamel.

Contraindication:

- Patients with respiratory problems, as the dust mist may aggravate their clinical condition.
- Patients who have kidney failure problems or who undergo hemodialysis.
- Hypertensive and cardiac patients, for these cases it will be advisable to seek medical advice before applying the product.
- Suspend use of the product in case of allergic reactions.

Instructions for Use:

- Apply a lip balm for greater patient comfort, avoiding dryness and cracks.
- Start the prophylaxis procedure following the manual of the blasting equipment manufacturer, and according to this procedure, fill the reservoir with sodium bicarbonate powder up to the inner step of the reservoir.
- Keep a distance of 3 to 5 mm from the tip of the equipment and the tooth surface, always with the angle towards the tooth, in order not to injure the gingival tissues.
- Special care should be taken with patients who have restorations, as in some cases, polishing of the restored surface will have to be performed in order to reduce roughness. And with patients who have exposed dentin, for example, in cervical lesions.

Form of presentation

- Pack with 250 g.
- Flavors: Mint and strawberry.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As a chemical, bicarbonate should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

Only trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.

BICARBONATO DE SÓDIO
SODIUM BICARBONATE
BICARBONATO DE SODIO

REV.01

- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product after the expiration date or if the packaging is damaged or tampered with.
- Keep out of reach of children.
- Avoid contact with eyes. If there is contact, wash immediately with plenty of water;
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- Keep the product in its packaging always closed.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- Avoid contact with moisture.
- Store in a ventilated and dry place, away from direct sunlight and at room temperature.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Disposal:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Reuse:

The bicarbonate powder is single-use. Any amount remaining in the equipment reservoir after the procedure must be discarded and must not be returned to the original packaging.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

- Sodium bicarbonate (NaHCO_3)
- Molecular formula: NaHCO_3
- Molar mass: 84.01 g/mol
- Average particle size: 38 μm

Symbology:

Symbol	Symbol Title	Description
--------	--------------	-------------

BICARBONATO DE SÓDIO
SODIUM BICARBONATE
BICARBONATO DE SODIO

REV.01

	Validity	Indicates the date after which the health product should not be used.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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