ENGLISH**Indication:**

Disposable Micro Applicators are indicated for the controlled application of substances in dental procedures. They can be used for the deposition of disclosing solutions (caries, bacterial plaque, among others), adhesive agents, sealants, eugenol, hemostatic solutions, fluorides and materials used in teeth whitening techniques, even helping in the movement of the gel during the procedure. The bendable shaft facilitates access to different regions of the oral cavity and favors precision in application.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing dental procedures that involve topical application of solutions, agents or fluid materials, regardless of age or clinical condition, as indicated and judged by the professional in charge.

Principle of operation:

Its function is based on the controlled transfer of the material by means of a flexible tip with fibers that retain small amounts of solution by surface tension, without significant absorption, allowing precise dosing and avoiding waste.

Contraindications:

- Do not use in patients with known hypersensitivity to the components of the device.
- Do not use with highly caustic substances or strong solvents that may compromise the integrity of the fibers or the stem.
- Do not use for the application of high viscosity or very thick materials that may deform or damage the fibers of the microapplicator.

Instructions for use:

- Remove an Applicator from the packaging.
- Touch the active tip to the solution to be applied.
- Apply the solution in the region of interest.
- The professional can choose the type of applicator (medium/blue, fine/green and extra fine/lilac), according to the size of the region to be worked and the amount of product to be applied.

Form of presentation:

- **Non-sterile product.**
- Single use.
- Package containing 100 units.
- Available in three sizes: medium, fine and extra fine.

BRUSH - APLICADORES DESCARTÁVEIS
BRUSH - DISPOSABLE MICRO APPLICATORS
BRUSH - MICROAPLICADORES DESECHABLES

REV.01

Cleaning procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure:

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and precautions:

- Product for professional use only.
- Keep out of reach of children.
- Do not subject the device to high temperatures.
- During operation, the professional must use all personal protective equipment, including goggles.
- Use eye protection on the patient.
- LOT number: keep the lot number shown on the product label; it must be provided to the manufacturer when reporting any complaint or incident.

Conservation and storage:

- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

Disposal:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

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REV.01

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Polypropylene and polyamide fibers.

Symbols:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates a medical device intended for one use, or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the medical device as defined in Regulation (EU) 2017/745 (MDR).
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	Prescription only	Sale and use restricted to qualified professionals or subject to prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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