CALBRADOR DE FOLGA INTERPROXIMAL INTERPROXIMAL GAP GAUGE CALBRADOR DE ESPACIO INTERPROXIMAL

REV.01

ENGLISH

Indication:

Interproximal space calibrators are indispensable tools to obtain accurate measurement of interproximal spaces resulting from interdental wear – stripping – made with sanding strip holders, diamond discs or any other abrasive device, in orthodontics. In addition, it corrects the dental discrepancy, calculated in the cephalometric study of orthodontic cases, and also acts in the manufacture of prostheses on implants in the measurements of prosthetic spaces.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing orthodontic treatment who require interproximal wear (stripping) to correct dental discrepancies, dental alignment, or optimize the available space in the arch. It also includes patients who are candidates for prosthetic rehabilitation on implants, in whom accurate measurement of prosthetic spaces is essential for treatment planning and execution.

Fundamentals of its operation and action:

Kit with six stainless steel blades of different thicknesses:

0.05; 0.10; 0.20; 0.30; 0.40 and 0.50 millimeters, providing precision and safety for the professional and comfort for the patient.


Contraindication:

- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and should not be used.
- Do not use cold sterilization methods (chemical process) for the sterilization of the Calibrator. These chemical agents can cause corrosion of the products.

Instructions for use:

Insert the Blade with the desired measurement to measure the size of the interproximal spaces resulting from interdental wear.

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Manual Cleaning Procedure:

- Use an area specifically reserved for cleaning.
- Rinse the calibrator under running water and keep them immersed in a bath of enzymatic detergent solution, following the manufacturer's instructions.
- After soaking and keeping it immersed, brush well away from the body with the neutral cleaning agent/enzymatic detergent. Care should be taken to avoid spreading contaminants from the spray during brushing procedures. Do not use wire brushes.
- Special care should be taken to clean cracks and other hard-to-reach areas. Visually inspect to confirm debris removal. Repeat the cycle if necessary.
- Rinse the calibrator thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towel.

Ultrasonic Cleaning Procedure

- Prepare an enzymatic detergent solution, following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place the calibrator in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use
- Completely immerse the device in the enzyme cleaner/detergent solution and program the equipment for at least fifteen (15) minutes.
- Rinse the calibrator thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towel.

Sterilization Requirements:

It is recommended to sterilize the product in a pressure-saturated steam autoclave at 121°C for 20 minutes, properly packaged in its own autoclave packaging, or the sterilization parameters should be followed in accordance with the relevant policies/rules of the country where the device is being used. It is the professional's sole responsibility to keep the autoclave equipment up to date to ensure the correct sterilization of the product.

Warnings and Precautions:

Only properly trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.

- Before use, inspect if the product is in perfect condition of use;
- Non-sterile product.
- Sterilize the product before use.
- During operation, the professional must use all personal protective equipment, including goggles.

CALBRADOR DE FOLGA INTERPROXIMAL
INTERPROXIMAL GAP GAUGE
CALBRADOR DE ESPACIO INTERPROXIMAL

REV.01

- Wear protective glasses on the patient.
- The Product is considered a contaminated instrument and proper handling precautions are required. Wear PPE, including gloves, masks, safety glasses, and a long-sleeved lab coat when handling, as is done with all other dental instruments.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- Keep the packaging in a dry place, away from light and away from the action of bad weather. Keep the products in their original packaging until the moment of use.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

The shelf life of the product may vary according to the conditions of use, frequency of sterilization, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

Additional information:

- 3R is not responsible if the instructions for use are not followed correctly, generating undesired results in the efficacy and safety of the product/patient.
- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

AISI 304 Stainless Steel

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.

CALIBRADOR DE FOLGA INTERPROXIMAL
INTERPROXIMAL GAP GAUGE
CALIBRADOR DE ESPACIO INTERPROXIMAL

REV.01

	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Medical device	Indicates that it is a medical device.
	Date of manufacture	Indicates the date the health product was manufactured.
	Maker	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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