ENGLISH

Indication:

Prophylaxis strip is indicated for clinical use in the removal of bacterial plaque and calculus of the interproximal region.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who require dental prophylaxis procedures, especially those with plaque accumulation or calculus in the interproximal region. It includes adolescents, adults and the elderly, as long as the procedure is indicated by the professional, aiming at promoting, maintaining or recovering periodontal health.

Fundamentals of its operation and action:

The prophylaxis strip is essential in prophylaxis procedures, as it easily accesses the proximal region. Its active face acts effectively in the removal of dental calculus.

The circular microblades, located around the holes on one of the sides of the strip, have as their main characteristic an action against calculus, preserving the integrity of the enamel and the restorative materials.

Contraindications:

- Marked tooth mobility.
- Do not use in areas with severe inflammation;
- Do not use the damaged or deformed product.
- Do not use on teeth at risk of fracture or fragile restorations.

Instructions for use:

Insert the strip into the interproximal space, making sure that the active face is facing the desired proximal face.

Hold the strip by the ends, pull it against the tooth and apply it in longitudinal movements in a back-and-forth form until the desired effect is obtained. Remove the strip.

Form of presentation

- Non-sterile product;
- One-time use
- Pack containing 10 units.
- Models: Micro-Blades and Soft Removal
- 4 mm

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

CINTA PROFILÁTICA
PROPHYLAXIS STRIP
CINTA PROFILACTICA

REV.01

- Professional use.
- Its reuse can cause cross-contamination.
- When there is a proven allergy to its component, the product should not be used.
- Contact of the strip with the gum may cause injury to the patient.
- Any serious incident related to the device must be reported to the manufacturer and the competent authority of the Member State in which the users and/or patients are established.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer.
- Packages must not be exposed to the weather.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Stainless steel

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Medical device	Indicates that it is a medical device.

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CINTA PROFILÁTICA
PROPHYLAXIS STRIP
CINTA PROFILACTICA

REV.01

Symbol	Symbol Title	Description
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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