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ENGLISH

Indication:

Diamond discs are used in a prosthesis laboratory or in a dental office for the preparation and finishing of prosthetic parts made of acrylic, porcelain or metal. The manufacturer is not responsible for the use of the product outside the established indication for use.

Intended user:

Dental professionals, including dental surgeons and dental technicians, duly trained to handle rotating instruments used in the preparation and finishing of prosthetic parts.

Contraindication:

- It should not be used intraorally.
- Do not use the product above the recommended RPM.
- Cleaning agents that contain chlorine or chlorides as an active ingredient are corrosive to stainless steel and should not be used.
- Do not use cold sterilization methods (chemical processes), as these agents can corrode the product.

Instructions for Use:

- Unlock the low-speed micromotor fitting or on the handpiece;
- Fix the disc in the chuck;
- Drive the low-speed motor;
- Start using.

Select the Diamond Disc that will be used and check the RPM according to the table below:

Diameter ISO	Speed Recommended	Speed Maximum	Pressure Work
(1/10 mm)	RPM (min ⁻¹)	RPM (min ⁻¹)	Weight (g)
005-014	70.000 – 90.000	250.000	30 – 50
015-023	60.000 – 80.000	150.000	30 – 50
024-040	40.000 – 50.000	80.000	50 – 100
041-100	20.000 – 40.000	50.000	50 – 100
110-220	7.000 – 12.000	20.000	50 - 100

Form of presentation:

- Non-sterile product;
- Single-use product;
- Models: Double-sided and single-sided;

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REV.01

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

- Before use, inspect if the product is in perfect condition for use;
- During operation, the professional must use all personal protective equipment, including goggles.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- Keep the packaging in a dry place, away from light and away from the action of bad weather. Keep the products in their original packaging until the moment of use.

Storage Conditions

- Storage should be in dry and clean conditions at room temperature.
- Keep the product stored in its original packaging until the moment of use.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service Life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- 3R is not responsible if the instructions for use are not followed correctly, generating undesired results in the efficacy and safety of the product/patient.
- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Diamond and Stainless Steel.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.

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REV.01

Symbol	Symbol Title	Description
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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