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**ENGLISH**

**Indication:**

The goat hair brush is a rotating instrument, which provides greater shine on the restored composite resin or ceramic surface. It should be used in the final stage of finishing and polishing installed in a contra-angle with a low-speed micromotor, additionally it can be used with diamond paste (or aluminum oxide paste) to obtain even more satisfactory results. The referred brush already has a rod (CA) for fitting in the contra-angle, that is, it does not require the use of an accessory chuck. The activation of this instrument can be done through compressed air or electric motor.

**Intended user:**

Duly qualified dental professionals.

**Target population:**

The product is intended for use in dental patients undergoing restorative procedures that involve the finishing and polishing of composite resin or ceramic restorations, according to the clinical indication of the dental surgeon.

**Fundamentals of its operation and action:**

Polisher. The goat hair brush has its action in the polishing and finishing stage, it has a versatile format for application on all faces of the restored surface.

**Contraindications:**

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than composite resin or ceramics.
- Do not use the product for purposes other than those for which it was developed.

**Instructions for use:**

The goat hair brush is a rotating instrument, which provides greater smoothness and shine of resins and ceramics used by dental professionals installed in a contra-angle with a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor.

- Carefully open the package making sure that the piece is in perfect condition for use.
- Then attach the polisher rod to the set against angle + micromotor of low rotation, making sure it is well locked.
- Start the polishing procedure by applying a small amount of polishing paste to the surface of the resin.
- It is recommended that refrigeration be used between intermittent use and should not apply pressure to the tooth surface.
- If necessary, repeat the procedure more than once.

**Form of presentation**

- Non-sterile product;

**Manual Cleaning Procedure**

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.

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REV.01

- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

**Ultrasonic Cleaning Procedure**

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

**Sterilization Requirements**

- Sterilization by pressure saturated steam autoclave at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used. It is the professional's full responsibility to keep the autoclave equipment in ideal conditions of use to ensure the correct sterilization of the product.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

**Warnings and Precautions:**

- Only legally qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.
- Before use, inspect whether the goat hair brush is in perfect condition for use;
- Non-sterile product.
- Sterilize the product before use.
- Make sure the goat hair brush is properly seated in the low-speed motor.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Ensure that the patient wears protective goggles.
- LOT Number: It is important to keep the lot number, identified on the product label, so that it can be reported to the manufacturer in the event of any occurrence.
- The disposal of the material must follow the relevant local legislation.
- Do not use in patients who have allergic reactions to nickel.

**Conservation and storage:**

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

**Discard:**

- The disposal of the material must follow the relevant local legislation.

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REV.01

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

**Service life:**

**5 uses (5 reuses); however, it is essential that a check be carried out before each use to confirm that the goat hair brush is in perfect condition for use.**

**Additional information:**

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: [www.microdont.com.br](http://www.microdont.com.br), the most up-to-date version of the instructions for use corresponding to the product.

**Composition:**

Active Part – Non-Abrasive Natural Bristles

Stem – Nickel-plated brass

**Symbology:**

| Symbol                                                                              | Symbol Title                           | Description                                                                                                      |
|-------------------------------------------------------------------------------------|----------------------------------------|------------------------------------------------------------------------------------------------------------------|
|  | Expiry date                            | Indicates the date after which the medical device should not be used.                                            |
|  | Non-sterile                            | Indicates a health product that has not undergone a sterilization process.                                       |
|  | Keep dry                               | Indicates a health product that needs to be protected against moisture.                                          |
|  | Medical device                         | Indicates that it is a medical device.                                                                           |
|  | Recyclable packaging                   | The packaging is recyclable.                                                                                     |
|  | See instructions for use               | Indicates the user's need to refer to the instructions for use.                                                  |
|  | Date of manufacture                    | Indicates the date the health product was manufactured.                                                          |
|  | Keep out of sunlight                   | Indicates a health product that needs protection from light sources.                                             |
|  | Manufacturer                           | Indicates the manufacturer of the medical device.                                                                |
|  | Lot Code                               | Indicates the manufacturer's batch code so that the batch or lot can be identified.                              |
|  | Catalog Number                         | Indicates the manufacturer's reference number so that the health product can be identified.                      |
| Rx Only                                                                             | Prescription only                      | Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation. |
|  | Unique Device Identification (UDI)     | Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).                              |
|  | QR Code – Access to the electronic IFU | QR Code to access the electronic version of the product's Instructions for Use (eIFU).                           |

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