ENGLISH

Indication:

The polishing brush impregnated with silicon carbide or alumina, is indicated for the finishing and polishing of dental restorations, especially those made of composite resins and ceramic materials. As it contains the abrasive incorporated into the bristles, it dispenses with the use of polishing pastes, providing uniform, efficient and easy application polishing.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the finishing and polishing of restorations in composite resin or ceramic materials, according to the clinical indication of the dental surgeon.

Fundamentals of its operation and action:

The silicon carbide or alumina polishing brush is a rotating instrument, installed in a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor. The silicon carbide or alumina polishing brush has its action on all types of restorative materials, available in brush and cup formats.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than those provided for in the indication for use.
- Do not use the product for purposes other than those for which it was developed.

Instructions for Use:

- Choose the shape of your preference, for the tooth surface in question.
- Unlock the low-rpm micromotor fitting.
- Insert the rod into the fitting, lock the rod into the low-rpm micromotor.
- Check that the rod is properly seated in the low-speed micromotor and that the coolant is working properly.
- Bring the micromotor with the brush to the tooth surface.
- Drive the micromotor at low speed without applying pressure so as not to generate excessive heat and move slowly evenly, if necessary, repeat the procedure.

Form of presentation:

- Non-sterile product;
- Single-use product;

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure:

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

ESCOVA POLIMENTO - CARBETO DE SILÍCIO / ALUMINA
POLISHING BRUSH - SILICON CARBIDE / ALUMINA
CEPILLO DE PULIDO - CARBURO DE SILICIO / ALUMINA

REV.01

Only specialized professionals in the dental medical field should use this product, as they have the scientific information necessary for its proper use.

- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product if the packaging is damaged, opened or wet.
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Silicon carbide impregnated polishing brush: The product is composed of CLA Ø 3/16" Brass rod with silicon carbide

Chemical composition of Brass CLA

Copper (Cu) -----55%

Lead (Pb) -----2.5 %

Zinc (Zn) -----Remaining

Silicon carbide

Nylon-type resin (nylon bristles)

ESCOVA POLIMENTO - CARBETO DE SILÍCIO / ALUMINA
POLISHING BRUSH - SILICON CARBIDE / ALUMINA
CEPILLO DE PULIDO - CARBURO DE SILICIO / ALUMINA

REV.01

Alumina impregnated polishing brush: The product is composed of CLA Ø 3/16" Brass rod with aluminum oxide.

Chemical composition of Brass CLA

Copper (Cu) -----55%

Lead (Pb) -----2.5 %

Zinc (Zn) -----Remaining

Chemical composition of aluminum oxide:

Aluminum (Al) -----0.07% max.

Oxygen (O)-----2.00% max.

Nylon-type resin (nylon bristles)

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the lot or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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