LENÇOL DE BORRACHA
RUBBER DAM
GOMA DIQUE

REV.01

ENGLISH

Indication:

The dental rubber sheet is indicated for absolute isolation of the operative field in order to contain saliva, blood or undesirable fluids in the place where the dental procedure will be performed. Maintaining asepsis conditions, facilitating the visualization of the operative field and keeping the field clean and dry.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing dental procedures that require absolute isolation of the operative field, such as adhesive restorations, endodontic treatments, pin cementation, indirect restorative procedures, moisture control in cavity preparations, and other interventions that require a clean, dry, and fluid-free field. The device can be used in patients of different age groups, according to the evaluation and indication of the professional in charge.

Fundamentals of its operation and action:

The dental rubber sheet is used for absolute isolation in dental procedures, product made of natural latex has flexibility and resistance, ensuring safety and effectiveness. It separates the tooth from contact with the mucous membranes and the tongue, creating an ideal operative field for the execution of various procedures in dentistry.

Contraindication:

- Contains natural latex, which may cause allergy.
- Do not use in patients who have reported sensitivity to latex.

Instructions for Use:

- Select the clamp to be used and/or the tie indicated for the tooth;
- Demarcate and drill the rubber sheet in the position of the teeth to be isolated (make sure that the drill is adjusted with the correct size hole, this is compatible with the tooth to be isolated);
- Mount the rubber sheet to the arch;
- If necessary, lubricate the edges of the sheet hole;
- Install the assembly to the tooth with the help of a clamp holder.

Form of presentation

- Non-sterile product;
- Single use;
- Package containing 26 units;
- Available in two sizes: 12.7 cm × 12.7 cm and 15 cm × 15 cm.
- Scent options: neutral or mint.

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

LENÇOL DE BORRACHA
RUBBER DAM
GOMA DIQUE

REV.01

Warnings and Precautions:

Only properly trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.

- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product if the packaging is damaged, opened or wet.
- After use, it should be discarded avoiding reuse, even if it still contains unused undrilled areas, thus avoiding the risk of cross-contamination.
- Forbidden to reuse.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 5 years, and must be used within this period to ensure its effectiveness and safety.

Life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

The product is manufactured in mint and natural versions.

Mint version

Natural latex: 97%

Vulcanizing: 2%

Natural Pigment: 1%

Mint flavouring: 1%

LENÇOL DE BORRACHA
RUBBER DAM
GOMA DIQUE

REV.01

Natural version

Natural latex: 97%

Vulcanizing: 2%

Natural Pigment: 1%

Composition of natural pigment:

Overseas blue pigment 15 (silicic acid, aluminum sodium salt, sulfur)

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates a health product intended for use or for use on a single patient during a single procedure.
	Keep dry	It indicates a health product that needs to be protected against moisture.
	Medical device	It indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	Contains or has the presence of natural rubber latex	Indicates the presence of natural rubber or natural rubber latex as the manufacturing material of the health product or its packaging.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	Indicates that the medical device must be used exclusively by a qualified health care professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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