MICRO DAM BLUE – BARREIRA GENGIVAL FOTOPOLIMERIZÁVEL
MICRO DAM BLUE – LIGHT CURING GINGIVAL BARRIER
MICRO DAM BLUE – BARRERA GINGIVAL FOTOPOLIMERIZABLE

REV.01

ENGLISH

Indication:

The gingival barrier is indicated for the protection of gum and mucous tissues during dental procedures that involve the use of bleaching agents, conditioning acids, or other potentially irritating chemicals. It can also be used for local isolation in specific areas of the oral cavity, helping to control humidity and delimit work areas during clinical procedures.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing dental procedures that involve the use of bleaching agents, conditioning acids, or other chemicals that may cause irritation to soft tissues.

Fundamentals of its operation and action:

The gingival barrier acts by forming a protective layer on the gingival and mucosal tissues adjacent to the teeth. After its application, the material is light-cured, forming a rigid and adherent layer, which works as a physical barrier against direct contact with chemical agents used in dental procedures, such as whitening gels or acids. This layer prevents the diffusion of irritating substances, reduces the risk of chemical burns and maintains the integrity of the soft tissues, ensuring safety and comfort during treatment.

Contraindication:

- Suspend use of the product in case of allergic reactions.
- Do not use on ulcerated areas;

Instructions for use:

- Position the lip, tongue and cheek retractor;
- Dry the gums and teeth corresponding to the region you want to protect;
- Apply a thin layer of the product (about 1 mm thick) on the gums and gingival papillae in such a way that you have total coverage and a good field of protection for the procedure you wish to perform. If diastemas are present, fill the space by applying the barrier through the tongue, in order to protect the interdental papillae and lingual gingival tissue;
- Polymerize the gum barrier for 20 to 30 seconds with a light curing. It is recommended to proceed with polymerization soon after the application of the resin;
- Check that the product covers the entire surface of the gums and that it is completely polymerized. To observe the upper part of papillae and gums, use a clinical mirror looking from incisal to cervical. Avoid testing it by forcing or trying to lift it because this can cause infiltration of the irritant;
- To remove the product, simply detach it forward with the help of an explorer. This should come off easily in its entirety or in large pieces.

Form of presentation:

- Low viscosity;
- Syringes with 2 g.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

MICRO DAM BLUE – BARREIRA GENGIVAL FOTOPOLIMERIZÁVEL
MICRO DAM BLUE – LIGHT CURING GINGIVAL BARRIER
MICRO DAM BLUE – BARRERA GINGIVAL FOTOPOLIMERIZABLE

REV.01

Sterilization Procedure:

Not applicable. As it is a chemical product, the gingival barrier should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.

- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product after the expiration date or if the packaging is violated.
- Avoid applying the product to the gums of patients who have recently undergone any type of oral surgical procedure that has any type of gum irritation.
- Keep out of reach of children.
- Avoid contact with eyes. If there is contact, wash immediately with plenty of water;
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- Keep the product in its packaging always closed.
- Store the product at temperatures from 5°C to 28°C
- Avoid contact with moisture
- Transport temperature 60°C for 5 days
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Shelf life:

Single-use product, with a shelf life restricted to the moment of clinical application. After being removed from the original packaging and used, any remaining amount of the product should be discarded immediately.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

MICRO DAM BLUE – BARREIRA GENGIVAL FOTOPOLIMERIZÁVEL
MICRO DAM BLUE – LIGHT CURING GINGIVAL BARRIER
MICRO DAM BLUE – BARRERA GINGIVAL FOTOPOLIMERIZABLE

REV.01

Composition:

GUM BARRIER

FUNCTION	COMPONENT	CAS RATING	PERCENTAGE (BY MASS)
Monomers	Urethane Dimethacrylate	72869-86-4	40 a 50%
	Hydroxyethyl Methacrylate	868-77-9	10 a 20%
Inert Charge	Nanometric Silica	69012-64-2	20 a 30%
Pigments	Patent Blue	129-17-9	1 a 3%
Photoinitiator	Canforaquinone	10373-78-1	0,5 a 1,0%
Co-initiator	Tertiary Amine: 4-ethyl dimethylaminobenzoate	10287-53-3	0,5- 2,0%

TIP

Stem: Stainless Steel

Barrel: Polypropylene

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Maker	Indicates the manufacturer of the health product as defined in Directive 90/385EED, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

MICRO DAM BLUE – BARREIRA GINGIVAL FOTOPOLIMERIZÁVEL
MICRO DAM BLUE – LIGHT CURING GINGIVAL BARRIER
MICRO DAM BLUE – BARRERA GINGIVAL FOTOPOLIMERIZABLE

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