## ENGLISH

**Indication:**

The product is indicated for the temporary sealing of cavities during restorative or endodontic treatments that have not been completed, as well as in situations requiring provisional closure on an emergency basis.

**Intended user:**

Duly qualified dental professionals.

**Target population:**

Intended for patients requiring provisional sealing of cavities during restorative or endodontic procedures that have been temporarily interrupted, including emergency care.

**Principles of operation and action:**

On contact with the patient's saliva, the paste hardens and seals the cavity, preventing contamination of the site by the oral environment (oral fluids and microorganisms).

**Contraindications:**

- Do not use in patients with hypersensitivity to any component of its formulation.
- Do not use as a base material under definitive restorative materials.

**Instructions for use:**

- Dry and isolate the application site;
- Apply the amount of filling material required to seal the cavity;
- Remove any excess so as not to interfere with occlusion;
- Remove the isolation and allow the product to come into contact with the oral environment;
- The setting time is 15 to 30 minutes after contact with saliva.

**Presentation:**

- Product supplied as a ready-to-use paste.
- Package containing 25 g.

**Cleaning procedure:**

Not applicable. This is a chemical product and therefore must not undergo any type of cleaning procedure.

**Sterilization procedure:**

Not applicable. As it is a chemical product, the temporary filling material must not be sterilized or subjected to any decontamination method.

**Warnings and precautions:**

- Only dental surgeons should use this product, as they have the scientific information required for its proper use.
- Wear PPE, including gloves, masks, safety glasses and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Use protective eyewear on the patient as well.
- Do not use the product after its expiry date or if the package has been tampered with.
- Keep out of the reach of children.
- Avoid contact with the eyes. In case of contact, rinse immediately with plenty of water.
- BATCH NUMBER: it is important to keep the batch number, shown in the product title, so that it can be reported to the manufacturer in the event of any occurrence.

**Conservation and storage:**

- Keep the product in its package, always closed.
- Store the product at a temperature between 10 °C and 30 °C.

**MICRO FILL - OBTURADOR PROVISÓRIO**  
**MICRO FILL - TEMPORARY FILLING MATERIAL**  
**MICRO FILL - OBTURACIONES PROVISÓRIAS**

REV.01

- Avoid contact with moisture.
- The product must be stored in a ventilated, dry place protected from light, so as to maintain ideal storage and transport conditions as well as its physical and chemical integrity.
- Keep the product in its original package until the moment of use. Inspect the package before use.
- The carrier must be informed of the contents. The care to be taken during transport is shown on the package, thereby ensuring protection of the product from dispatch until delivery to the customer. Packages must not be exposed to the weather.

**Disposal:**

- Disposal of the material must comply with the applicable local legislation.

**Expiry date:**

The product has a shelf life of 2 years and must be used within this period to ensure its efficacy and safety.

**Service life:**

Single-use product, with a service life restricted to the moment of clinical application. Once removed from the original package and used, any remaining quantity must be discarded immediately.

**Other additional information:**

- If necessary, 3R provides the instructions for use in printed format, at no additional cost, upon request by e-mail (qualidade@3r.ind.br) or by telephone: (+55 11-5525-0590 / +55 0800-773-5525), with a response within up to 7 business days of the request.
- Always check the website [www.microdont.com.br](http://www.microdont.com.br) for the most up-to-date version of the instructions for use corresponding to the product.

**Composition:**

- Zinc Oxide
- Calcium Sulfate
- Non-toxic Fluid Silicone

**Symbols:**

| Symbol                                                                              | Symbol Title                 | Description                                                                                                       |
|-------------------------------------------------------------------------------------|------------------------------|-------------------------------------------------------------------------------------------------------------------|
|  | Expiry date                  | Indicates the date after which the medical device must not be used.                                               |
|  | Keep dry                     | Indicates a medical device that needs to be protected from moisture.                                              |
|  | Medical device               | Indicates that the item is a medical device.                                                                      |
|  | Consult instructions for use | Indicates the need for the user to consult the instructions for use.                                              |
|  | Date of manufacture          | Indicates the date on which the medical device was manufactured.                                                  |
|  | Manufacturer                 | Indicates the manufacturer of the medical device.                                                                 |
|  | Batch code                   | Indicates the manufacturer's batch code so that the batch or lot can be identified.                               |
|  | Catalogue number             | Indicates the manufacturer's catalogue number so that the medical device can be identified.                       |
|  | For professional use only    | Indicates that the medical device is to be used exclusively by a qualified healthcare professional.               |
| Rx Only                                                                             | Sale by prescription         | Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation. |

**MICRO FILL - OBTURADOR PROVISÓRIO**  
**MICRO FILL - TEMPORARY FILLING MATERIAL**  
**MICRO FILL - OBTURACIONES PROVISÓRIAS**

REV.01

| Symbol                                                                            | Symbol Title                           | Description                                                                            |
|-----------------------------------------------------------------------------------|----------------------------------------|----------------------------------------------------------------------------------------|
|  | QR Code – Access to the electronic IFU | QR Code to access the electronic version of the product's Instructions for Use (eIFU). |

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