ENGLISH

Indication:

It is indicated for cementation and fixation of metal prosthetic parts such as crowns and bridges, in addition to acting in the cementation of bands in orthodontic appliances. Fluidity is its main characteristic, which enables the proper placement of prosthetic and orthodontic components.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who require cementation of metal prosthetic parts — such as crowns, bridges, or pins — or who use orthodontic bands as part of their treatment. Any patient requiring fixed prosthetic rehabilitation or orthodontic treatment that involves banding is included.

Fundamentals of its operation and action:

GLASS IONOMER C is a radiopaque glass ionomer cement, chemically activated used for the cementation of crowns, bridges, inlays, onlays, orthodontic appliances. It chemically adheres to metal and tooth structure. It is a cement whose liquid contains polycarboxylic acids that guarantee the ability of chemical adhesion to the tooth structure. It has biocompatibility, releases fluoride for anticariogenic action and has excellent adhesion to the tooth structure.

Contraindication:

- The product is contraindicated for patients with allergies to any of the components of the formula.
- Do not use on a damp surface.
- Do not use in cementing laminated veneers and gluing brackets.

Instructions for Use:

- Thoroughly clean the enamel, dentin, and metal surfaces of prosthetic parts, drying them thoroughly without dehydrating the tooth.
- Make sure that there are no eugenol residues, as they interfere with the adhesion of the cement; It is recommended to use eugenol-free temporary cements.
- Smooth the internal surfaces of the restoration with a finishing carbide or fine-grained diamond drill bit.
- Place the powder and liquid on the glass plate or paper pad for mixing.
- Add the powder to the liquid in small portions, gradually increasing until you get the right consistency.
- When dripping the liquid, always keep the bottle in an upright position.
- Mixing time: 0 to 1 minute.
- Setting time: 3 to 4 minutes.
- Optimum handling temperature: 18 °C to 32 °C (above this range decreases working time; below it increases).
- In deep cavities, the use of protective varnish is recommended.
- The cement will be at the correct consistency when it is possible to apply it to the preparation and the inner surface of the restoration with a rigid brush.
- Avoid contact with water or saliva during insertion and setting.
- Apply an even layer of approximately 1 mm.
- Position the restoration or band and remove excess before solidification is complete.
- Immediately rinse metal instruments with cold water to prevent adhesion and hardening of cement residue.

MICROGLASS - IONÔMERO DE VIDRO C
MICROGLASS - GLASS IONOMER C
MICROGLASS - IONÓMERO DE VIDRIO C

REV.02

Work Accessories:

- Glass plate or paper pad, plastic spatula and measuring spoon.

Ratio for mixing:

- 1 level measure of powder to 2 drops of liquid

Form of presentation

- Powder: 10g
- Liquid: 8ml

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As it is a chemical, the glass ionomer should not be sterilized nor subjected to any method of decontamination.

Warnings and Precautions:

- Avoid contact of the material with water and saliva during insertion and setting.
- Eugenol cements should not be used in temporary restorations.
- Phenolic substances such as eugenol affect the adhesion of the material.
- Do not use in cementing laminated veneers and bonding brackets.
- Do not use after the expiration date.
- Keep out of reach of children.
- For dental use only.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

Single-use product, with a service life restricted to the moment of clinical application. After being removed from the original packaging and used, any amount remaining in the device should be discarded immediately.

Storage and conservation:

- Store the product in a dry place, at temperatures between 15 and 30 °C. After use, keep the jar closed.
- Store away from medicines such as: Phenol, Eugenol, Paramonochlorophenol etc.
- Transport the package tightly closed in a dry and clean place. Do not exceed a temperature of 37°C and a relative humidity of 75%.

Additional information:

MICROGLASS - IONÔMERO DE VIDRO C
MICROGLASS - GLASS IONOMER C
MICROGLASS - IONÓMERO DE VIDRIO C

REV.02

- If necessary, 3R provides instructions for use in printed form, at no additional cost, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

PHASE 1: POWDER	
RAW MATERIAL	FUNCTION
Fluorine Strontium Aluminosilicate Glass	Restorative Agent

PHASE 2: LIQUID	
RAW MATERIAL	FUNCTION
Polyacrylic acid copolymer	Binding agent

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or subject to prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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