ENGLISH

Indication:

It is indicated for lining cavities and sealing cracks.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures who require cavity lining.

Fundamentals of its operation and action:

The GLASS IONOMER F adheres to the enamel and dentin. It can be applied to dentin without any protection and still be treated with the acid etching technique in the same way that tooth enamel is treated. It works as a great base material for composites and amalgams. As a rule, no special retention is required due to chemical adhesion to the cleaned enamel and dentin. It is radiopaque.

Contraindication:

- The product is contraindicated for patients with allergies to any of the components of the formula.
- Do not use on a damp surface.

Instructions for Use:

- Thoroughly clean the enamel, dentin, and metal surfaces of prosthetic parts, drying them thoroughly without causing dental dehydration.
- Make sure that there are no eugenol residues, as they impair the adhesion of the cement; The use of eugenol-free temporary cements is recommended.
- Smooth the internal surfaces of the restoration with a finishing carbide drill or fine-grained diamond drill bit.
- In very deep regions, with a risk of pulp exposure, a cement based on calcium hydroxide can be applied beforehand.
- In deep cavities, apply a layer of protective cement and, after setting, proceed with the definitive restoration in composite resin or amalgam.
- Before handling, gently shake the powder bottle to avoid compaction. Deposit the powder (removing the excess) and liquid on the glass plate or pad of paper.
- Add the powder to the liquid in small portions, gradually incorporating until it reaches a homogeneous consistency suitable for application.
- When dripping the liquid, always keep the bottle in an upright position; The first drop can be discarded due to the possible presence of air bubbles.
- Mixing time: 0 to 1 minute.
- Setting time: 3 to 4 minutes.
- Optimum handling temperature: 18 °C to 32 °C (above this range, the working time decreases; below, it increases).
- The correct consistency is achieved when the cement can be applied to the dental preparation and the inner surface of the restoration using a rigid brush.
- Avoid any contact with water or saliva during insertion and the setting period.
- Apply an even layer of approximately 1 mm.
- Immediately rinse metal instruments with cold water to prevent adhesion and hardening of cement residue.

Work Accessories:

MICROGLASS - IONÔMERO DE VIDRO F
MICROGLASS - GLASS IONOMER F
MICROGLASS - IONÓMERO DE VIDRIO F

REV.01

Glass plate or paper pad, plastic spatula and measuring spoon.

Ratio for mixing:

- 1 level measure of powder to 1 drop of liquid

Form of presentation

- Powder: 10g
- Liquid: 8ml

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As it is a chemical, the glass ionomer should not be sterilized nor subjected to any method of decontamination.

Warnings and Precautions:

- Avoid contact of the material with water and saliva during insertion and setting.
- It is recommended to use calcium hydroxide-based products if the cavity to be restored from the remaining dentin is less than 1.5 mm.
- Eugenol cements should not be used in temporary restorations.
- Phenolic substances such as eugenol affect the adhesion of the material.
- Temporary restorations should be made with eugenol-free cements.
- During insertion and setting, avoid contact of water and saliva with the material.
- Do not use after the expiration date.
- Keep out of reach of children.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- For dental use only

Discard:

The disposal of the material must follow the relevant local legislation.

Expiration date:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Shelf life:

After opening, the powder and the liquid may be used until the expiration date, provided the bottles are kept tightly closed and stored under the conditions indicated below. Any unused mixed portion must be discarded immediately.

Conservation and storage:

- Store the product in a dry place, at temperatures between 15 and 30 °C. After use, keep the jar closed.
- Store away from medicines such as: Phenol, Eugenol, Paramonochlorophenol etc.
- Transport the package tightly closed in a dry and clean place. Do not exceed a temperature of 37°C and a relative humidity of 75%.

MICROGLASS - IONÔMERO DE VIDRO F
MICROGLASS - GLASS IONOMER F
MICROGLASS - IONÓMERO DE VIDRIO F

REV.01

Other additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional cost, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product

Composition:

Liquid: Polyacrylic Acid

Powder: Calcium Barium Fluorosilicate Aluminum.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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