ENGLISH**Indication:**

The GLASS IONOMER R is indicated for restorations in deciduous teeth in all cavitory classes (I, II, III, IV and V), and can be used as both provisional and definitive restorative material. It is also indicated for temporary restorations in permanent teeth in the same classes. In addition, it can be used as a basis for restorations with different materials and for the sealing of pits and fissures. The product performs well in the treatment of non-carious lesions, such as abfractions, and in the Atraumatic Restorative Treatment (ART) technique.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for pediatric and adult patients who require temporary or permanent restorations in primary teeth and temporary restorations in permanent teeth. It also includes patients with non-carious lesions, such as abfractions, and those indicated for treatment using the ART technique.

Fundamentals of its operation and action:

The GLASS IONOMER R is a self-curing cement based on glass ionomer, developed for fast-setting dental restorations, presented in two phases: powder and liquid. It has high biocompatibility, releases fluoride with anticariogenic action and provides excellent chemical adhesion to the tooth structure.

Contraindication:

- The product is contraindicated for patients with allergies to any of the components of the formula.
- Do not use on a damp surface.

Instructions for Use:

- After cavity preparation, perform prophylaxis with pumice stone and water. Use absolute isolation; if this is not possible, perform effective relative isolation to avoid saliva contamination.
- Condition the cavity with polyacrylic acid (using the liquid of the R GLASS IONOMER) for 15 seconds. Then, rinse thoroughly and dry carefully, avoiding dehydration of the tooth substrate.
- Smooth the inside surfaces of the cavity with a finishing carbide bur or a fine-grit diamond bur. It is not necessary to perform mechanical retentions, as the cement has chemical adhesion to the tooth.
- Shake the powder bottle to avoid compaction. Place the powder (removing the excess) and liquid on the paper mixing pad, adding the powder to the liquid gradually, starting with small portions and increasing as the mixture develops.
- When dripping the liquid, keep the bottle in an upright position to ensure accuracy; It is recommended to discard the first drop due to the possible presence of air bubbles.
- Mix for 30 seconds to 1 minute. The optimum handling temperature is 18 °C to 32 °C (temperatures above reduce working time; temperatures below prolong it).
- Introduce the cement into the cavity using a syringe or spatula, avoiding the incorporation of air bubbles.
- The setting time of the material is 3 to 4 minutes; avoid contact with saliva or moisture during this time.

Working Accessories:

- Glass plate or paper pad, plastic spatula and measuring spoon.

Mixing ratio:

MICROGLASS R - IONÔMERO DE VIDRO - A2 / A3
MICROGLASS R - GLASS IONOMER - A2 / A3
MICROGLASS R - IONÓMERO DE VIDRIO - A2 / A3

REV.01

- 1 level measure of powder to 1 drop of liquid

Form of presentation

- Powder: 10g
- Liquid: 8ml

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As it is a chemical, the glass ionomer should not be sterilized nor subjected to any method of decontamination.

Warnings and Precautions:

- Avoid contact of the material with water and saliva during insertion and setting.
- It is recommended to use calcium hydroxide-based products if the cavity to be restored from the remaining dentin is less than 1.5 mm.
- Eugenol cements should not be used in temporary restorations.
- Phenolic substances such as eugenol affect the adhesion of the material.
- Temporary restorations should be made with eugenol-free cements.
- During insertion and setting, avoid contact of water and saliva with the material.
- Do not use after the expiration date.
- Keep out of reach of children
- LOT Number: It is important to keep the lot number, identified on the product label, to be reported to the manufacturer in the event of any occurrence.
- For dental use only

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

Single-use product, with a shelf life restricted to the moment of clinical application. After being removed from the original packaging and used, any amount remaining in the device should be discarded immediately.

Conservation and storage:

- Store the product in a dry place, at temperatures between 15 and 30 °C. After use, keep the jar closed.
- Store away from medicines such as: Phenol, Eugenol, Paramonochlorophenol etc.
- Transport the package tightly closed in a dry and clean place. Do not exceed a temperature of 37°C and a relative humidity of 75%.

Additional information:

MICROGLASS R - IONÔMERO DE VIDRO - A2 / A3
MICROGLASS R - GLASS IONOMER - A2 / A3
MICROGLASS R - IONÓMERO DE VIDRIO - A2 / A3

REV.01

- If necessary, 3R provides instructions for use in printed form, at no additional cost, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Powder: Sodium fluoroaluminosilicate and phosphate salts, iron oxide when staining is required.

Liquid: Polyacrylic copolymer solution.

Symbology:

Symbol	Symbol Title	Description
	Expiry date	Indicates the date after which the health product should not be used.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or subject to prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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