ENGLISH

Indication:

Surgical grade paper is indicated for sterilization and maintenance of sterility. It is commonly used in dental-medical-hospital environments for the packaging of instruments and/or articles that will be submitted to sterilization processes, constituting one of the links in the chain of protection of professionals and patients.

Intended user:

Qualified professionals in the dental, medical and hospital areas, including dental surgeons, doctors, nurses and technicians/assistants responsible for the preparation, packaging and sterilization of materials.

Fundamentals of its operation and action:

The surgical grade paper is presented in the form of a reel or envelope, consisting of a disposable package composed of a thermo-resistant plastic sheet and another of SMS (spunbond-meltblown-spunbond) both fused laterally, where there are indicators that show when it has been subjected to heat and autoclave pressure. The indicators will change color after sterilization to indicate that the envelopes have been processed, this proof when used in conjunction with routine biological indicators (spore testing), provides a greater level of sterility assurance.

Instructions for Use:

- To present the envelope, proceed as follows:

- Carefully place the items in the envelope, if it is sharp, we recommend that you protect this region with gauze;
- Remove the liner from the adhesive strip;
- Fold the adhesive paper flap;
- Press the paper flap firmly and seal the envelope securely.

- For the coil presentation, proceed as follows:

- Cut a piece of the reel according to the size of the item you will sterilize, so that the packages remain longitudinal and never transverse, because the longest sealing edges will be the original factory ones, that is, the edges less sensitive to breaking/opening.
- With the sealer already heated and ready for use, stretch the package well so as not to form wrinkles and seal one end;
- Carefully place the items in the package, if it is sharp, we recommend that you protect this region with gauze;
- Seal the other end, leaving about 3 cm beyond the sealed region so the package can be opened aseptically;

Contraindications:

Reusing the packaging is prohibited. Items should not be squeezed in the package or envelope.

Warnings and Precautions:

- Before use, inspect whether the envelope is suitable for use.
- Never use if the envelope is damaged.
- Single-use product
- The indicator does not guarantee sterilization, but it is essential that it is present in each package, as this will avoid the use of a package that has not been autoclaved;
- Keep out of reach of children;
- Do not use the product after the expiration date;
- Product for professional use only.

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REV.01

- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Conservation and storage:

Store the envelopes in the original packaging, closed in a clean, cool place, away from light and away from the action of bad weather.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 5 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Paper: Surgical grade paper 60g/m2 (ISO11607-1-2 and EN868-3)

Film: Transparent film, light blue/green, PET/PP, 52µm (EN868-5)

Sterile Indicator: Imported ink; Chemical sterile indicator printing size ≥100mm2, no fading or color loss after sterile.

Symbology:

Symbol	Symbol Title	Description
	Use-by date	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Attention	Do not use the device if it is damaged.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.

This review cancels and replaces the previous revision

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REV.01

Symbol	Symbol Title	Description
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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