DENTAL PROPHYLACTIC PASTE
PASTA PROFILÁCTICA DENTAL
PASTA PROFILÁTICA ODONTOLÓGICA

REV.01

ENGLISH

Indication:

Prophylactic paste is indicated in dental prophylaxis for the removal of extrinsic stains, plaque, materia alba and oral debris.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing dental prophylaxis as clinically indicated by the dental surgeon.

Fundamentals of its operation and action:

The prophylactic paste has an adequate consistency to be placed on the teeth without running into the oral cavity. The product can be applied either with Robinson's brush or with rubber cups for prophylaxis. Can be used on adults and children. It is a product that is abrasive enough to remove all types of residues from the tooth surface without causing undue abrasion of the enamel, dentin or cementum. In addition to acting as a cleaning agent, it gives the hard tissue an aesthetic and highly polished appearance.

Contraindication:

- This product should not be used in case of hypersensitivity (allergy) to any of the components of the formula.

Instructions for Use:

- Dispense the necessary amount of paste into a dappen jar.
- Then press the brush/bowl onto the paste.
- Apply a small amount to the surface of the teeth, using low rotation until complete removal of stains and plaque. Then rinse.
- Avoid keeping the brush/cup on the tooth surface without the paste, as this, combined with continuous (non-intermittent) use of the low-speed handpiece, results in excessive heating of the tooth surface.

Form of presentation

- 90 g pack.
- Contains fluoride.
- Available in mint or tutti-frutti flavors.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As a chemical product, the prophylactic paste should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

Only properly trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.

- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product after the expiration date or if the packaging is violated.

DENTAL PROPHYLACTIC PASTE
PASTA PROFILÁCTICA DENTAL
PASTA PROFILÁTICA ODONTOLÓGICA

REV.01

- During the procedure, care must be taken to avoid excessive abrasion of any restorative material present, these are susceptible to wear due to their low hardness.
- Apply carefully to areas with exposed dentin.
- Do not use the product continuously or for prolonged periods, as this will cause the tooth enamel to wear down.
- Keep out of reach of children.
- Do not ingest. If ingested, it may cause nausea.
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- Keep the product in its packaging always closed.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

- The disposal of the material must follow the relevant local legislation.

Expiry date:

The product has an expiry date of 2 years from the date of manufacture, and must be used within this period to ensure its effectiveness and safety.

Shelf life:

The paste is intended for multi-patient use within the expiry date stated on the packaging. The portion dispensed into the dappen jar is single-use: any remaining amount after the procedure must be discarded immediately and must not be returned to the original packaging.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

The product is manufactured in the Tutti-Frutti and Mint versions.

Tutti-Frutti flavour version:

Substance Quantity %

Carboxymethylcellulose ----- 1.30

Sorbitol 70% ----- 14,00

DENTAL PROPHYLACTIC PASTE
PASTA PROFILÁCTICA DENTAL
PASTA PROFILÁTICA ODONTOLÓGICA

REV.01

Sodium Lauryl Sulfate-----	5.50
Fumed silica -----	1.50
Calcium Carbonate-----	25.00
Quartz 325-----	16.00
Methyl Paraben-----	0.15
Sodium Saccharin -----	0.30
Aroma de Tutti-Frutti-----	0,40
Deionized water -----	31.15
Bi-distilled Vegetable Glycerin -----	4.60
Sodium Fluoride-----	0.10
	100,00

Mint flavour version:

Substance Quantity %

Carboxymethylcellulose -----	1.30
Sorbitol 70% -----	14,00
Sodium Lauryl Sulfate-----	5.50
Fumed silica -----	1.50
Calcium Carbonate-----	25.00
Quartz 325-----	16.00
Methyl Paraben-----	0.15
Sodium Saccharin -----	0.30
Mint flavouring-----	0.40
Deionized water -----	31.15
Bi-distilled Vegetable Glycerin -----	4.60
Sodium Fluoride-----	0.10
	100,00

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.

DENTAL PROPHYLACTIC PASTE
PASTA PROFILÁCTICA DENTAL
PASTA PROFILÁTICA ODONTOLÓGICA

REV.01

Symbol	Symbol Title	Description
	Catalogue number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorised representative in the European Community.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or subject to prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code - Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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