
ENGLISH

Indication:

The diamond spiral polisher is a rotating instrument, which helps in the finishing and polishing stage of the restored surface in composite resin. It allows greater versatility, since its configuration can act on all sides of the dental surface. It is a two-step system: Initially, the pink polisher (medium grain) should be used for wear, roughness removal and other initial corrections and then the gray polisher (fine grain) to obtain greater shine and smoothness. It must be installed in a contra-angle with a low-speed micromotor, and the activation of this can be done through compressed air or electric motor.

Pay attention to the green ring present in the polisher that identifies it as a composite resin polisher according to the image below:

Intended user:

Duly qualified dental professionals.

Target population:

Patients who need composite resin restorations and who require adequate finishing and polishing of these surfaces to ensure aesthetics, smoothness, comfort and longevity of the restorative treatment.

Fundamentals of its operation and action:

Abrasion. The spiral polisher has its action in the polishing and finishing stage.

Contraindications:

- It should not be used on restorative materials other than composite resin.
- Do not use the polisher counterclockwise.
- Do not operate at a speed higher than 15,000 RPM.

Instructions for use:

- Carefully open the package making sure that the piece is in perfect condition for use.
- Then attach the polisher rod to the low rotation counter angle + micromotor assembly, making sure it is well locked. It is of paramount importance that the direction of rotation is set to rotate clockwise (also identified with the letter "F" on the micromotor).
- The recommended rotation range per minute is from 8,000 to 12,000 RPM, and you should not exceed 15,000 RPM under any circumstances.
- Start polishing the restored surface with the pink spiral (medium grain)

SPINGLOSS - POLIDORES ESPIRAIS PARA RESINA
SPINGLOSS - SPIRAL POLISHERS FOR RESIN
SPINGLOSS - PULIDORAS EN ESPIRAL PARA RESINA

REV.01

- After finishing use with the pink spiral, wash the tooth surface thoroughly, so that no residue of the polisher previously used is present at the time it is replaced by gray (fine graining).
- It is recommended that refrigeration be used between intermittent use and should not apply pressure to the restored tooth surface.
- If necessary, repeat the procedure more than once.

Form of presentation

- Non-sterile product;
- Kit with 2 polishers
- Pink: Medium; Grey: Fine

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by pressure saturated steam autoclave at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used. It is the professional's full responsibility to keep the autoclave equipment in ideal conditions of use to ensure the correct sterilization of the product.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Warnings and Precautions:

Only legally qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.

SPINGLOSS - POLIDORES ESPIRAIS PARA RESINA
SPINGLOSS - SPIRAL POLISHERS FOR RESIN
SPINGLOSS - PULIDORAS EN ESPIRAL PARA RESINA

REV.01

- Do not exceed the limit of 5 uses and always check beforehand that the polisher is in optimal condition.
- Before use, inspect whether the spiral polisher is in perfect condition for use;
- Non-sterile product.
- Sterilize the product before use.
- Make sure that the spiral polisher is properly seated in the low-speed motor.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Ensure the patient wears protective goggles.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

5 uses (5 reuses), however, it is of paramount importance that before each use a check is carried out, that is, if the dental polishers are in perfect condition for use

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Metal rod composed of stainless steel.

Rubber polisher with abrasive composed of diamond particles.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.

SPINGLOSS - POLIDORES ESPIRAIS PARA RESINA
SPINGLOSS - SPIRAL POLISHERS FOR RESIN
SPINGLOSS - PULIDORAS EN ESPIRAL PARA RESINA

REV.01

Symbol	Symbol Title	Description
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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