POLIGLOSS E SPEEDYGLOSS
POLIGLOSS AND SPEEDYGLOSS
POLIGLOSS Y SPEEDYGLOSS

REV.01

ENGLISH

Indication:

The polishing parts (cup, torpedo and disc) are rotating instruments, which provide greater smoothness and shine of resins, used by dental professionals, and are installed in a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the finishing and polishing of composite resin restorations, as indicated by the dentist.

Fundamentals of its operation and action:

Abrasion. The product is composed of silicone tips mounted on a metal rod, acts by polishing and finishing, and has varied shapes for application on different tooth surfaces.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than the composite resin.
- Do not use the product for purposes other than those for which it was developed.

Instructions for use:

- Start polishing the resin on the most abrasive piece and change it until you reach the least abrasive piece (in the following order: Coarse, Medium, Fine);
- Fit the piece to the metal rod.
- Unlock the low speed micromotor fitting;
- Insert the rod into the socket.
- Lock the rod on the low rpm micromotor;
- Check that the rod is properly fitted to the low speed micromotor and that the coolant is working properly;
- Apply polishing paste to the polishing part;
- Take the micromotor with the piece for polishing to the tooth surface;
- Drive the micromotor at low speed without applying pressure so as not to generate excessive heat and move slowly evenly;

Form of presentation:

- Non-sterile product;
- Models: Cup, torpedo and disc
- Kit with 18 polishers + 1 chuck;

KIT SPEEDYGLOSS:

POLIGLOSS E SPEEDYGLOSS

POLIGLOSS AND SPEEDYGLOSS

POLIGLOSS Y SPEEDYGLOSS

REV.01

KIT POLIGLOSS:

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by pressure saturated steam autoclave at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Warnings and Precautions:

POLIGLOSS E SPEEDYGLOSS

POLIGLOSS AND SPEEDYGLOSS

POLIGLOSS Y SPEEDYGLOSS

REV.01

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- Before use, inspect whether the cup, torpedo and rubber disc are in perfect condition for use;
- Non-sterile product.
- Sterilize the product before use.
- Make sure the rubber bowl/torpedo/disc is properly seated in the low-speed motor.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Place protective goggles on the patient.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

5 uses (5 reuses), however, it is of paramount importance that before each use a check is carried out, that is, if the dental polishers are in perfect condition for use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Silicone

Abrasive: Aluminum Oxide

Snap-on mandrel

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.

POLIGLOSS E SPEEDYGLOSS
POLIGLOSS AND SPEEDYGLOSS
POLIGLOSS Y SPEEDYGLOSS

REV.01

	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's lot code so that the lot or batch can be identified.
	Catalogue number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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