ENGLISH

Indication:

Ultrasonic Burs are used in conjunction with dental ultrasound equipment and are indicated for various situations, including calculus removal, cavity preparation and finishing, post decementation and endodontic procedures. The different models allow the professional to choose the most appropriate for each clinical situation.

- **Scaler Tips** – In general, they are indicated for the removal of supragingival calculus.
- **Periodontics Tips** – They are indicated for the removal of supragingival and/or subgingival calculus on the entire surface of the tooth and several other clinical situations detailed below, in the section corresponding to periodontics.
- **Cavitory Preparation Tips** - In general, they are indicated for the removal of decayed tissue without damaging the adjacent teeth, the specific indications are detailed below, in the section corresponding to the tips for preparation.
- **Endodontics Tips** - In general, they are indicated for the removal of calcifications, pulp nodules, fractured instruments and root canal cleaning, the specific indications for each situation are described in the section corresponding to endodontics.

• **A/B/C**

Function: Periodontal Treatment/Prophylaxis

Application: The 3 tips are indicated for action in supragingival areas, assisting in the removal of clinically visible calculus on the buccal and lingual/palatine surfaces and if applicable in interproximals. Each tip has a different angle to better adapt to the clinical situation.

• **G1 / G2 / GD1 / GD2**

Function: Periodontal Treatment/Prophylaxis

Application: **G1:** It is indicated to assist in the removal of clinically visible stones present in the crown limiting to the gingival margin in small amounts and/or little adherent. It has a tip with a thinner active part, and for this reason it can be used on the proximal surfaces, however, if appropriate to the clinical situation, it can also be used on the buccal and lingual/palatine surfaces.

G2: It is indicated to assist in the removal of clinically visible bulky stones, present in the crown and limited to the gingival margin. Ideally applicable on buccal and lingual/palatine surfaces, as its tip with a straight active part is suitable for these indications. For the tip to be effective, it is recommended that it be in contact with the surface of the tooth where it is necessary to act.

• **G3 / G4 / GD3 / GD4**

Function: Periodontal Treatment/ Prophylaxis

Application: **G3 and G4:** They are indicated to assist in acting in supragingival interproximal and subgingival areas (shallow pockets), as it has an active tip with a pointed shape.

• **G5/ GD5**

Function: Periodontal Treatment/Prophylaxis

Application: It has a tip with a thinner active part and with a hole for high cooling, indicated to act in the removal of clinically visible supragingival stones, poorly adhered and/or in small quantity, ideally indicated on the proximal surfaces, however if appropriate to the clinical situation it can also be used on the buccal and lingual/palatal surfaces.

• **G6/ GD6**

Function: Periodontal Treatment/Prophylaxis

Application: It is indicated to assist in the removal of clinically visible bulky stones ideally present on buccal and lingual/palatine surfaces, as its straight active tip is suitable for these indications. For the tip to be effective, it is recommended that it be in contact with the tooth surface.

- **P1 / P3 / PD1 / PD3**

Function: Periodontal Treatment

Application: **P1:** It has a finer and slightly sharp active tip, it is indicated to assist in the removal of subgingival stones, with apical formation at the gingival margin, usually difficult to visualize clinically.

P3: It has a thinner and more elongated active tip, it is indicated to assist in the removal of subgingival stones, with apical formation at the gingival margin, usually difficult to visualize clinically.

- **P3D / PD3D**

Function: Periodontal Treatment

Application: It has an active diamond tip indicated to assist in the treatment of the root surface/regularization after periodontic surgeries and removal of sharp edges of restored crowns.

- **P4/PD4**

Function: Periodontal and Endodontic Treatment

Application: It has an active tip with a tapered end and with high cooling, it helps in the removal of calculus in interdental regions, and if appropriate to the clinical situation it can also be used on other faces. In addition to the indication in periodontics, this tip can be used in endodontics to act on the pulp chamber, refining the access surgery and removing the roof of the pulp chamber.

- **P4D / PD4D**

Function: Periodontal and Endodontic Treatment

Application: Indicated to assist in periodontal procedures where a tip with tapered end, diamond and with high cooling is necessary, in addition, it can be used in endodontics to act in the pulp chamber, refining the access surgery and removing the roof of the pulp chamber.

- **E1 / ED1**

Function: Endodontic Treatment

Application: It is indicated in the support for files used during cleaning and widening of channels, with an angle of 120°, commonly used in anterior teeth.

- **E2 / ED2**

Function: Endodontic Treatment

Application: It is indicated in the support for files used during cleaning and widening of channels, with an angle of 90°, commonly used in posterior teeth.

- **E3 / ED3**

Function: Endodontic Treatment

Application: It is indicated to assist in endodontic procedures for condensation and gutta-percha cutting, besides facilitating the removal of pins, breaking the cementation line between the pin and the canal walls.

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REV.01

Function: Endodontic Treatment

Application: It is indicated to reach the pulp chamber and assist in the removal of the pulp chamber roof, refinement of access surgery, internal calcifications and pulp nodules.

- **E5 / ED5**

Function: Endodontic Treatment

Application: Indicated for cutting and condensing gutta percha.

- **E5D / ED5D**

Function: Endodontic Treatment

Application: Indicated to assist in the removal of calcifications, particularly useful for removal while smoothing the walls, it also assists to act in isthmus areas.

- **ME1**

Function: Endodontic Treatment

Application: Indicated to perform PUI (Passive Ultrasonic Irrigation), and assist in work in isthmus or hard-to-reach areas.

- **ME2D**

Function: Endodontic Treatment

Application: Indicated to assist in the refinement of access surgery, removal of pulp chamber roof, removal of nodules, calcifications or when a lateral action is necessary.

- **ME3D**

Function: Endodontic Treatment

Application: Indicated to assist in the location of calcified canals, refinement of access surgery, removal of the pulp chamber ceiling, removal of nodules and calcifications.

- **ME4**

Function: Endodontic Treatment

Application: Indicated to act in isthmus areas and assist in the removal of fractured files.

- **ME5**

Function: Endodontic Treatment

Application: Indicated to assist in the location of calcified channels, unfilling and removal of fractured files.

- **ME6D**

Function: Endodontic Treatment

Application: Indicated to assist in the removal of the pulp chamber roof, refinement of access surgery, location of canals, removal of coronary obstructions, removal of restorative materials, removal of calcifications, removal of temporary and permanent cements and removal of fiber pins.

- **ME7D**

Function: Endodontic Treatment

Application: Indicated to assist in the refinement of access surgery, removal of the pulp chamber ceiling, removal of nodules and calcifications.

- **ME8**

Function: Endodontic Treatment

Application: Indicated to assist in the beginning of unfilling.

- **ME10**

Function: Endodontic Treatment

Application: Indicated for cutting and condensing gutta percha.

- **ME12**

Function: Endodontic Treatment

Application: Indicated to assist in the removal of metal pin.

- **ME18**

Function: Endodontic Treatment

Application: Indicated to assist in the location of calcified channels, removal of fractured file and act in isthmus areas.

- **MR1**

Function: Endodontic Treatment

Application: Indicated to assist in retreatments, in the defilling stage.

- **SB1 / SBD1**

Function: Cavitory Preparation

Application: Active tip with diamond sphere that assists in access and adjustments of occlusal and cervical surfaces, in addition to the removal of hard materials in restoration procedures.

- **SB2 / SBD2**

Function: Cavitory Preparation

Application: Tip with active diamond part (half of the upper sphere) to assist in the preparation and adjustments of proximal surfaces (mesial) without damaging the adjacent tooth.

- **SB3 / SBD3**

Function: Cavitory Preparation

Application: Tip with active diamond part (half of the lower sphere) to assist in the preparation and adjustments of proximal surfaces (distal) without damaging the adjacent tooth.

- **SBL / SBDL**

Function: Cavitory Preparation

Application: Active diamond part (left half sphere), with an angle of 45° to the left, allows access in difficult interproximal areas, assisting in the removal of caries lesions without damaging adjacent areas.

- **SBR/SBDR**

Function: Cavitory Preparation

Application: Active diamond part (right half sphere), with an angle of 45° to the right, allows access in difficult interproximal areas, assisting in the removal of caries lesions without damaging adjacent areas.

• C1 / C2 / CD1 / CD2

Function: Uncementing Pins

Application: **C1:** Assists in the removal of metal pins. Used to transfer ultrasonic power to cements whose retention is done by micromechanical principles. Not recommended for use in adhesive type cements. It does not have irrigation.

C2: Assists in the removal of metal pins. Used to transfer ultrasonic power to cements whose retention is done by micromechanical principles. Not recommended for use in adhesive type cements. It has irrigation.

• C3 / CD3

Function: Uncementing Pins

Application: It is indicated to use it on the cementation line, with the aim of assisting in the decementation of metal pins. It does not have a hole for irrigation.

• Equipment Compatibility

G1, G2, G3, G4, G5, G6, P1, P3, P3D, P4, P4D, E1, E2, E3, E3D, E5, E5D, SB1, SB2, SB3, SBL, SBR, C1, C2, C3, ME1, ME2D, ME3D, ME4, ME5, ME6D, ME7D, ME8, ME10, ME12, ME18, MR1

Adiel Plus | Advance SE | DMC | DTE | Easy Sonic | GG Dent | Gnatus | J Morita | Olsen | Ortus
Import | Satelec | Schuster | Spartan | Ultradent | VDW | VT150

GD1, GD2, GD3, GD4, GD5, GD6, PD1, PD3, PD3D, PD4, PD4D, ED1, ED2, ED3, ED3D, ED5, ED5D, SBD1, SBD2, SBD3, SBDL, SBDR, CD1, CD2, CD3

ALT | Advance S | Emissionic (MMO) | Kavo | Kavo Unik Komfort | Kondortec | Mectron | Microdont
| Mini Piezon | Ortus | PM 100 | PM 200 | WH | Woodpecker

Fundamentals of its operation and action:

It assists in various procedures such as supra and subgingival scaling, de-cementation of pins, removal of decayed tissue and endodontic procedures.

Instructions for Use:

- The surgeon is responsible for checking the products before use to know if they are suitable for their intended purpose.
- Select the correct maximum power as indicated by the ultrasound equipment.
- Choose tip/insert
- Insert the tip with the help of the wrench in the ultrasound, check if the equipment has the correct cooling and apply it in the place according to your indication.

Manual Cleaning Procedure:

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning;
- Wash with the help of a soft-bristled brush and enzymatic detergent under running water;
- Submerge the product in a specific cleaning solution for surgical instruments (follow the manufacturer's instructions);
- Remove from solution;
- Care should be taken to clean cracks and other hard-to-reach areas. Visually inspect to confirm debris removal. Repeat the cycle if necessary.
- Rinse the ends thoroughly under running water until they are visibly clean.
- Pat the ends dry with clean compressed air or paper towels.

Ultrasonic Cleaning Procedure:

For ultrasonic cleaning, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution/enzymatic detergent following the agent manufacturer's instructions for correct concentration, exposure time, temperature, and water quality.
- Place the dental tips in the appropriate location for the ultrasonic equipment, following the manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse the dental tips thoroughly under running water until they are visibly clean.
- Pat the ends dry with clean compressed air or paper towels.

Sterilization Requirements:

This item applies to tips purchased in non-sterile form and to the reuse of products purchased sterile.

- Pack in material suitable for sterilization;
- Sterilize in an autoclave at 121 °C for 20 minutes, or follow the sterilization parameters according to the applicable policies and regulations in the country of use.

Contraindications:

- Cleaning agents with chlorine or chloride as an active ingredient are corrosive to stainless steel and should not be used.

Warnings and Precautions:

- Non-sterile product
- Sterilize the product before use.
- Sterilization is not replaced by cleaning.
- Once sterilized, the package should be opened only at the time of use, under aseptic conditions.
- Exclusive Product for dental use
- Do not use the product if the packaging is damaged, opened or wet.
- Those with the suffix "D" at the end of the description are diamond products, that is, they produce wear and tear of dental materials and tissues. They have 5 reuses, and it is still necessary to check beforehand if it is in conditions of use.
- During the operation, the professional must use all personal protective equipment, including protective goggles. Wear protective goggles on the patient.
- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- 3R is not responsible for use of the product outside of the intended use.

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REV.01

- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- The disposal of the material must follow the relevant local legislation.

Conservation and storage:

- The product must be in an airy, dry place and away from light in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the products in their original packaging until the moment of use.
- The carrier must be informed of the contents.
- The precautions to be taken in transport are demonstrated in the packaging, thus ensuring protection of the instruments from shipment to delivery to the customer.
- Packaging should not be left to the elements.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

- 5 uses (5 reuses); however, it is essential that a check be carried out before each use to confirm that the ultrasonic tips are in perfect condition for use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional cost, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

AISI 440C Stainless Steel

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product.

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REV.01

Symbol	Symbol Title	Description
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	Indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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