ENGLISH

Indication:

The diamond bur is a device for use in dental treatment that is connected to high-speed instruments (high-speed handpiece) and eventually at low speed (micromotor coupled in contra-angle with adapter mandrel). It is used in several clinical situations, the most common being the performance of cavitary preparations, in addition to being used in the wear of a wide variety of materials found in dentistry. Dental materials are: restorative materials (hard metals, amalgam, porcelain, resin composites, glass ionomers, in addition to dental tissues, such as enamel).

3R Indústria E Comércio Ltda is not responsible when the product is used outside the intended use proposed.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in patients undergoing dental procedures that involve the controlled removal of dental structures through wear, according to the clinical need determined by the dental surgeon.

Applications of use by type of dental diamond tips:

The Dental Diamond Tips can be used in high-speed instruments (High- speed handpiece), always with water/air cooling to reduce the heat generated by friction, as well as to eliminate residues of dental tissue and/or materials.

The indication for the use and function of the device is related to the design and size of the active tip that will be used according to the clinical situations. Grain size and/or grain type are important specifications in diamond tips that ensure the interchangeability and compliance of these products with the requirements of dental practice.

Dental diamond tips have different granulations: extra coarse (181 µm), coarse (151 µm), medium (91 - 126 µm), fine (46 µm) and extrafine (30 µm). Extra coarse (black) and coarse (green) tips are used to quickly reduce the thick layers of the tooth structure. They promote a faster reduction than medium tips, which, in turn, are used in routine tasks involving enamel, dentin, and cavitary preparations in general. Fine tips (red color) and extra-fine tips (yellow color) are used in wall finishing and at the edges of cavity preparation in the surface finishing of direct or indirect composite resin or ceramic restorations

Dental diamond tips come in different shapes and sizes, and each is coded (ISO code) for identification.

				
806	31 + 4	001	524	006
Material of Working part (Diamond)	Shank and overall length	Form (round, cylinder, etc.)	Type (Toothing, coarse, etc.)	Nominal size Largest working part diameter in 1/10mm

	(FG / RA)			
ISO 6360-1	ISO 1797-1 / ISO 6360-1	ISO 6360-2	ISO 6360-3	ISO 2157
COD. ISO:	806.314.001.524.006			

Source: ISO 6360-1 Standard

The following dental diamond tips are examples of tips commonly used in cavity preparations:

- **Spherical:** Indicated for initial access to caries lesions. Other uses are recommended, according to the diameter of the device:
 - **Small Diameter:** minimally invasive cavity preparation and additional retentions in cavity preparation for silver amalgam restoration.
 - **Medium Diameter:** delimitation of orientation grooves in preparations for total crown, veneers and indirect restorations in general.
 - **Large diameter:** coronary opening, pulp chamber roof removal and beveling at the cavo-superficial angle of extensive cavitory preparations in direct and indirect restorations.
- **Spherical with collar:** Indicated for use in cavitory preparations, helping, for example, in the initial access to carious lesions, opening or retentions in direct dentistics, making orientation grooves in indirect restorations, and regularization of the cavo-superficial angle.
- **Conical:** Although the shape of these diamond tips is conical, they have variations that will condition the clinical application, described below:
 - **Inverted Conical:** Used primarily to create mechanical retention. They are useful in the preparation of self-retentive cavitory preparations in extensive restorations. In addition, larger ones flatten and regularize the gingival wall, pulp and dihedral angles.
 - **Conical torpedo:** Its tip crosses the proximal contact point in cavity preparations inside and outside the access. These tips are used on proximal faces for indirect metal, porcelain, or composite resin restorations.
 - **Conical Flat End:** Indicated to make cavitory preparations with expulsive walls and defined internal angles. It is also indicated for making niches for support of the removable partial prosthesis.
 - **Rounded end conical:** Indicated for making cavity preparation for porcelain/composite resin inlay and onlay with expulsive walls and rounded angles and onlay preparation with cusp protection ending in

conical chamfer. It is also indicated for the preparation of orientation grooves and reduction of the buccal face in prosthetic preparations of laminates.

- **Flat-ended cylindrical:** These are used to perform parallelism on the walls of the preparation and make the dihedral angles sharper. The fine tips are used in the wall finishing of the cavity preparation. Indicated for making and finishing the finish in the shape of a shoulder or marginal step in prosthetic preparations.
- **Wheels:** Indicated to create orientation grooves for wear of the buccal face.
- **Flame:** They help in the making of concave bevels, preparations on the lingual/palatine surfaces in anterior teeth, finishing in the cavo-superficial angle and can also be useful in wear for occlusal adjustment.

Contraindications:

- To avoid fracture and distortion of the instrument, do not apply excessive force when inserting dental diamond tips at high speed;
- Respect the working speed (rpm Max) defined on the product label, otherwise it may generate undesirable heat;
- Do not use dental diamond tips as levers and/or incorrect angles;
- Do not apply excessive pressure to the dental diamond tips, as this can cause undesirable heat and/or cause tip failures and risks to the patient;
- Chlorine or chloride cleaning agents as an active ingredient are corrosive to stainless steel and should not be used.
- Prolonged immersion of dental diamond tips in disinfectant solutions can result in corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) for sterilization of dental diamond tips. These chemical agents can cause corrosion of the products.

Instructions for use:

- Insert dental diamond tips only into a high-speed handpiece that meets the requirements of the applicable IEC 60601 series of standards and other relevant product standards. Apply the abrasive part to the surface to be prepared. Use dental diamond tips only for recommended indications.
- Select the dental diamond tips according to the speed (we have a speed table - RPM- for the use of the tips).

Diameter ISO	Recommended Speed	Maximum Speed	Working pressure
(1/10mm)	RPM (Rotation per Minute)	RPM (Rotation per Minute)	Weight (g)
008-018	250.000 – 300.000	450.000	20 – 50
021-031	100.000 – 120.000	160.000	20 – 50
033-047	30.000 – 90.000	120.000	20 – 50
049-063	10.000 – 25.000	80.000	20 – 50
065-090	5.000 – 15.000	45.000	20 - 50

Using:

PONTA DIAMANTADA
DIAMOND BURS
PUNTA DE DIAMANTE

REV. 01

- Move the dental diamond tips continuously to avoid localized heating and consequent breakage of the tips, causing risks to the patient;
- Use proper cooling methods to ensure proper cooling of the parts used for treatment;
- Use aspiration techniques that do not compromise the cooling of the parts involved in the treatment;
- Eye protection should be worn to protect against ejected particles.
- Used dental diamond tips should be considered contaminated instruments and appropriate handling precautions are required. Wear personal protective equipment (PPE), including gloves, masks, and safety glasses when handling, as with all other dental instruments

Form of presentation:

- Non-sterile product;
- Individual package containing 1 unit.

Cleaning Procedures:

- Clean and sterilize the dental diamond tips before each use as per cleaning and sterilization item
- It is recommended that cleaning be done immediately after the end of the procedure, to avoid drying out contaminants, which can make cleaning difficult.
- After cleaning, check the diamond tips to ensure that any traces of contamination have been removed. Repeat cleaning if necessary.
- The cleaning of dental diamond tips can be done in two ways: Manual or Ultrasonic.

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows.

- Use an area specifically reserved for cleaning;
- Wash the dental diamond tips under running water and keep them immersed in a bath of neutral pH cleaning solution (enzymatic detergent) following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent (enzymatic detergent). Care should be taken to avoid spreading contaminants by spraying during the brushing process. Do not use wire brushes.
- Care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal. Repeat the cycle if necessary.
- Wash the dental diamond tips thoroughly under running water until they are visibly clean.
- Dry the device with clean compressed air or paper towels.

Ultrasonic Cleaning Procedure

For ultrasonic cleaning, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent) following the agent manufacturer's instructions for correct concentration, exposure time, temperature, and water quality.
- Place the dental diamond tips in the appropriate location for the ultrasonic equipment, following the manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Wash the dental diamond tips thoroughly under running water until they are visibly clean.
- Dry the device with clean compressed air or paper towels.

Sterilization requirements for Dental Diamond Tips:

PONTA DIAMANTADA
DIAMOND BURS
PUNTA DE DIAMANTE

REV. 01

- It is recommended to sterilize the Dental Diamond Tips by saturated steam autoclave under pressure at 121°C for 20 minutes, or the sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used.
- It is the professional's sole responsibility to keep the autoclave equipment up to date to ensure the correct sterilization of the product.

USE LIMITATIONS: The recommended maximum number of reuses of the product is 5 use cycles (5 reprocessing cycles), however it is of paramount importance that before each use a check is carried out, that is, if the dental diamond tips are in perfect condition for use, otherwise they can overheat the dental tissues and lead to the pulp response triggering a painful sensation in the patient, in addition to not complying with the necessary wear requirement.

Warnings and Precautions:

- Clean and sterilize dental diamond tips before each use according to cleaning and sterilization items
- For user and patient safety, the dental diamond tip should only be used by qualified dental professionals.
- Before use, check that the dental diamond tips are in perfect condition for use, especially in relation to the diamond area and integrity of the rod
- Before use, make sure that the stem of the dental diamond tips is fully seated in the high-speed handpiece;
- Never force the dental diamond tips into the high-speed handpiece. In case of difficulty in fitting, check the turbine and the dental diamond tips.
- Respect the working speed (RPM Max) defined on the product label, otherwise you may generate undesirable heat.
- In cases of returning dental diamond tips to 3R (manufacturer), make sure that they have been effectively cleaned, sterilized and properly packaged for return.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed of any occurrences with the manufacturer.
- Material must be disposed of in accordance with local regulations.

Storage conditions:

Storage should be dry and clean at room temperature.

Keep the product stored in its original packaging until the moment of use.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Additional Information:

1. The reuse processes described above have been validated and are suitable and effective for the reprocessing of dental diamond tips. 3R is not responsible if the instructions for use are not properly followed, resulting in undesirable results in product/patient efficacy and safety.

2. If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request via email (qualidade@3r.ind.br) or telephone (11-5525-0590 / 0800-773-5525), and we will get back to you no later than 7 business days from the date of the request.

Always check the website: www.microdont.com.br, the most current version of the usage instruction corresponding to the product.

Composition:

PONTA DIAMANTADA
DIAMOND BURS
PUNTA DE DIAMANTE

REV. 01

AISI 303 stainless steel

Diamond

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Maker	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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				
806	31 + 4	001	524	006
Material of Working part (Diamond)	Shank and overall length (FG / RA)	Form (round, cylinder, etc.)	Type (Toothing, coarse, etc.)	Nominal size Largest working part diameter in 1/10mm
ISO 6360-1	ISO 1797-1 / ISO 6360-1	ISO 6360-2	ISO 6360-3	ISO 2157
COD. ISO:	806.314.001.524.006			

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