RÉGUA ENDODÔNTICA
ENDODONTIC RULER
REGLA ENDODÓNTICA

REV.01

ENGLISH

Indication:

Millimeter endodontic rulers are indicated for measuring the length of endodontic instruments (such as files, gutta-percha cones and absorbent paper cones), used during root canal treatment procedures.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing endodontic procedures that require accurate determination of the working length of the instruments used in root canal preparation and filling. It includes adults, adolescents and, when indicated by the professional, children who need endodontic treatment, always under the supervision of qualified professionals.

Fundamentals of its operation and action:

Endodontic rulers work by direct mechanical measurement, using a millimeter scale that allows adjusting and controlling the length of endodontic instruments and materials. This action ensures the accuracy of the working length, contributing to the safety and efficacy of endodontic treatment.

Instructions for use:

- Make sure the ruler is clean and dry before use.
- Position the endodontic instrument (file, gutta-percha cone, or absorbent paper) in the groove or guide of the ruler.
- Adjust the instrument's rubber cursor or silicone marker to the desired reference point.
- Read the length on the millimeter scale to determine the exact measurement of the instrument.
- Adjust the working length as planned for endodontic treatment.
- After use, clean and disinfect the ruler according to the manufacturer's recommendations.

Form of presentation:

- Non-sterile product
- Colors: Pink, Silver, Green, Red or Lilac

Manual Cleaning Procedure:

If manual cleaning is the only option available, proceed as follows.

- Use an area specifically reserved for cleaning;
- Rinse the endodontic ruler under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent) following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent (enzymatic detergent). Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Care should be taken to clean cracks and other hard-to-reach areas. Visually inspect to confirm debris removal. Repeat the cycle if necessary.
- Rinse the endodontic ruler thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towels.

Ultrasonic Cleaning Procedure:

For ultrasonic cleaning, it is recommended to proceed as follows:

RÉGUA ENDODÔNTICA
ENDODONTIC RULER
REGLA ENDODÓNTICA

REV.01

- Prepare a pH-neutral cleaning solution (enzymatic detergent) following the agent manufacturer's instructions for correct concentration, exposure time, temperature, and water quality.
- Place the endodontic ruler in the appropriate location for the ultrasonic equipment, following the manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse the endodontic ruler thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towels.

Sterilization requirements for the endodontic ruler:

- It is recommended to sterilize the endodontic ruler in an autoclave with saturated steam under pressure at 121°C for 20 minutes, or to follow the sterilization parameters according to the relevant policies/regulations of the country where the device is being used.
- It is the professional's sole responsibility to keep the autoclave equipment up to date to ensure the correct sterilization of the product.

Contraindications:

Not applicable

Warnings and Precautions:

- Before use, inspect if the product is in perfect condition of use;
- Non-sterile product.
- During operation, the professional must use all personal protective equipment.
- Ensure the patient wears protective eyewear.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- Keep the packaging in a dry place, away from light and away from the action of bad weather. Keep the products in their original packaging until the moment of use.

Conservation and storage:

- Store in a clean, dry place.
- Transport in a clean, dry place at room temperature.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

The service life of the product may vary according to the conditions of use, frequency of sterilization, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

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REGLA ENDODÓNTICA

REV.01

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website www.microdont.com.br for the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Aluminum 6351
Width: 12.7mm
Length: 80mm

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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