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ENGLISH

Indication:

Soft Bur is indicated for surgical procedures in soft tissues of the oral cavity and to assist in periodontal treatments, including gingival contouring (gingivectomy), exposure of intraosseous implants, removal of granulation tissue, treatment of gingival hyperplasias, papillectomy, exposure of impacted or semi-impacted teeth, exposure of deep cervical caries and removal of the gingival groove for impression.

Intended user:

Duly qualified dental professionals.

Target population:

It is intended for patients who require surgical interventions in oral soft tissues or periodontal and restorative procedures that involve gingival remodeling, removal of pathological tissues or exposure of dental structures.

Fundamentals of its operation and action:

Ceramics allow vessels to coagulate and reduce the tendency to hemorrhage.

It provides high stability, allowing it to act on the gums quickly and reliably. In addition, minimal warmth prevents necrosis during treatment.

Contraindication:

- Do not use Soft Bur on hard tissues, such as tooth enamel or bone tissue.
- Do not use the bur in conjunction with irrigation, as this compromises its proper performance.

Instructions for Use:

- Dry the operative region completely; The surface of the soft tissues should be clean and dry before use.
- Be sure to turn off the cooling of the high-speed pen. Do not use water spray during the procedure, as refrigeration can increase bleeding and impair cutting efficiency.
- Fit the tip to the high-speed pen as directed by the manufacturer, ensuring its proper locking.
- Activate the pen and apply the tip to the mucosa to be removed, using rotation between 300,000 and 500,000 rpm, preferably at maximum speed. Lower rotations can cause improper vibration, reduced efficiency, and risk of tissue damage.
- Insert the tip into the tissue and make the cut. In this rotation range, kinetic energy is converted into controlled heat at the moment of contact, promoting immediate coagulation. The tip does not heat up, as it is made of heat-resistant material; Heat is transferred directly to the tissue, producing instant hemostasis.

Form of presentation

- Non-sterile product;

Manual Cleaning Procedure:

- Use an area specifically reserved for cleaning.
- Rinse the tip under running water and keep it immersed in a bath of enzymatic detergent solution, following the manufacturer's instructions.
- After soaking and keeping it immersed, brush well away from the body with the neutral cleansing agent. Care should be taken to avoid spreading contaminants from the spray during brushing procedures. Do not use wire brushes.
- Special care should be taken to clean cracks and other hard-to-reach areas. Visually inspect to confirm debris removal. Repeat the cycle if necessary.
- Rinse the tip thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towel.

Ultrasonic Cleaning Procedure

- Prepare an enzymatic detergent solution, following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place the bur in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse the bur thoroughly under running water until it is visibly clean.
- Dry the device with clean compressed air or paper towel.

Sterilization Requirements:

It is recommended to sterilize the product in a pressure-saturated steam autoclave at 121°C for 20 minutes, properly packaged in its own autoclave packaging, or the sterilization parameters should be followed in accordance with the relevant policies/rules of the country where the device is being used. It is the professional's sole responsibility to keep the autoclave equipment up to date to ensure the correct sterilization of the product.

Precautions and Warnings:

- Before use, inspect that the tip is in perfect condition for use.
- Sterilize the product before use.
- It is important not to use the pen's refrigeration, as refrigeration can increase the tendency to hemorrhage, preventing blood clotting.
- Speeds below 300,000 rpm can cause injury
- During operation, the professional must use all personal protective equipment, including goggles.
- Wear protective glasses on the patient.
- Soft Bur is considered a contaminated instrument and proper handling precautions are required. Use PPE, including gloves, masks, and safety glasses when handling, as is done with all other dental instruments.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- Keep the packaging in a dry place, away from light and away from the action of bad weather. Keep the products in their original packaging until the moment of use.

Disposal:

- The disposal of the material must follow the relevant local legislation.

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REV.01

• Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

• 5 uses (5 reuses), however, it is of paramount importance that before each use the check is carried out, that is, if the zirconia bur is in perfect condition for use.

Conservation and storage:

• Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
 • The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

Additional information:

• If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
 • Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Ceramic and stainless steel tip

Symbology:

Symbol	Symbol Title	Description
	Use-by date	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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