TIRAS DIAMANTADAS DE AÇO GALVANIZADO
GALVANIZED STEEL DIAMOND STRIPS
TIRAS DIAMANTADAS DE ACERO GALVANIZADO

REV. 01

ENGLISH

Indication:

Stainless steel diamond strips are instruments used in dental procedures for finishing, polishing, and interproximal wear. They allow for controlled removal of excess restorative material and adjustment of the contact point between adjacent teeth.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative or orthodontic procedures (interproximal wear) as indicated by the dentist.

Applications:

Finishing and polishing of the proximal surfaces of restorations in composite resin and other dental materials, as well as for the removal of excesses, adjustment of contact points between adjacent teeth and interproximal wear.

Fundamentals of its operation and action:

Stainless steel diamond sanding strips act by abrasive mechanical action, promoting the controlled removal of excess restorative material or tooth enamel through abrasion between the active surface of the strip and the area to be worn. This action allows for interproximal wear, finishing, and polishing of restorations or dental surfaces.

Contraindications:

- Do not use in patients with known hypersensitivity to metals, especially stainless steel;
- Do not use in areas with severe inflammation or active bleeding;
- Do not use the product that is damaged, deformed or has lost its abrasive coating.

Instructions for use:

Stainless steel diamond sanding strips feature a neutral center (non-abrasive) that is used to insert the sandpaper into the proximal faces without wearing down the contact area:

- Hold the sandpaper with the pulps of your thumb and index fingers.
- Insert the sandpaper with the neutral center into the interproximals of the teeth.
- Move from the buccal to palate/lingual and vice versa, from the cervical region towards the occlusal surface and vice versa.
- To facilitate insertion on the proximal surface, a lubricant can be used with sandpaper.

Form of presentation

- Non-sterile product;
- One-time use
- Pack containing 6 units.

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

TIRAS DIAMANTADAS DE AÇO GALVANIZADO
GALVANIZED STEEL DIAMOND STRIPS
TIRAS DIAMANTADAS DE ACERO GALVANIZADO

REV. 01

Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- Goggles should be worn to protect against ejected particles.
- Wear PPE, including gloves, masks, safety goggles, and long-sleeved lab coats when handling the strip, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Stainless Steel and Diamond.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure

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REV. 01

	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community
	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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