BANDA MATRIZ DE AÇO INOX
STAINLESS STEEL MATRIX BAND
BANDA DE MATRIZ DE ACERO INOXIDABLE

REV. 01

ENGLISH

Indication:

Steel matrix bands are auxiliary devices used in the reconstruction of dental walls and in the protection of adjacent teeth during cavity preparation, avoiding wear beyond the element of interest. They are indicated to promote the proximal contour and reestablish the point of contact, presenting adequate rigidity so that they do not move during the condensation of the amalgam or the insertion of restorative materials, such as composite resin.

Intended user:

Duly qualified dental professionals.

Target population:

The steel matrix bands are intended exclusively for dental surgeons and other qualified dental professionals for use in restorative procedures.

Fundamentals of its operation and action:

They wrap around the tooth forming a wall.

Contraindications:

- Patients with known hypersensitivity to stainless steel or any component of the material.
- Do not use in areas with severe inflammation or active bleeding;
- Cases in which the tooth does not present minimum structural conditions for retention or stabilization of the matrix band.

Instructions for use:

Matrix bands replace one or more missing tooth walls and are ideal for promoting proximal contouring and maintaining the contact point, and are rigid enough not to shift during amalgam condensation or insertion of restorative material such as composite resin.

Stainless steel matrix band in roll:

Cut the matrix band to the desired size

Fit the cut-out piece into the "matrix retainer"

Wrap the tooth with the matrix band

Fit to desired size

For a better adaptation, use "wedges" in the interproximals of the tooth.

Perform the restoration

Remove the "wedges"

Release the "matrix retainer"

Remove the matrix band

Stainless steel matrix band in strips:

Fit the strip to the "matrix retainer"

Wrap the tooth with the matrix band

Fit to desired size

For a better adaptation, use "wedges" in the interproximals of the tooth

Perform the restoration

Remove the "wedges"

BANDA MATRIZ DE AÇO INOX
STAINLESS STEEL MATRIX BAND
BANDA DE MATRIZ DE ACERO INOXIDABLE

REV. 01

Release the "matrix retainer"
Cut the matrix band
Remove the matrix band

Tofflemire matrix band

Fit the Tofflemire into the "matrix retainer"
Wrap the tooth with the matrix band
Fit to desired size
For a better adaptation, use "wedges" in the interproximals of the tooth
Perform the restoration
Remove the "wedges"
Release the "matrix retainer"
Cut the matrix band
Remove the matrix band

Sectional matrices

- Carefully open the packaging of the matrix bands, making sure that all parts are inside the package.
- Perform absolute isolation of the operative field.
- Select the appropriate size of the matrix according to the clinical situation (occlusal-gingival height).
- Position the matrix in the cervical direction and accommodate it in the region that will be restored with the concave surface facing the tooth to be restored.
- If the matrix cannot be accommodated, pre-wedge or lightly wear the prepared tooth with a sanding strip to allow the matrix to fit without being dented.
- If pre-wedging has been performed, remove the wedge to allow the matrix to be seated, then reposition it as needed according to the space.
- With the specific clamp for fixing the matrix, fix it on the tooth surface. In extreme cases, the matrix conformation can be adjusted with the aid of a burnisher or exploration probe. However, be careful at this stage not to crease the matrix, which already has the ideal shape.
- Once the procedure is complete, remove the clamp, then the wedge and, finally, the matrix.

WITH EXTENSION

WITHOUT EXTENSION

Form of presentation

- Non-sterile product;
- One-time use

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

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REV. 01

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

- Before use, inspect that the Matrix Band is in perfect condition for use.
- During operation, the professional must use all personal protective equipment, including goggles.
- Wear protective glasses on the patient.
- The Matrix Band is considered a contaminated instrument and proper handling precautions are required. Wear PPE, including gloves, masks, safety glasses, and a long-sleeved lab coat when handling it, as with other dental instruments.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Storage conditions:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Stainless Steel

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the medical device should not be used.

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REV. 01

Symbol	Symbol Title	Description
	Medical device	Indicates that it is a medical device.
	CE Mark	Indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	For professional use only	Indicates that the medical device must be used exclusively by a qualified healthcare professional.
	Date of manufacture	Indicates the date the medical device was manufactured.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the medical device can be identified.
	Authorized representative in the European Community	Indicates the authorized representative in the European Community
	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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