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ENGLISH

Indication:

The Polyester Matrix Band are films that assist in the protection and reconstruction of tooth walls, they are ideal to protect the adjacent tooth during steps pertinent to acid conditioning and/or adhesive application, in addition they can assist in class IV restorations performed freehand, they also promote proximal contour and maintain the point of contact.

3R is not responsible for improper use of the product, especially when used for purposes other than those for which it was specifically developed.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures.

Fundamentals of its operation and action:

They surround the tooth, forming a wall that allows the light from the light-curing device to pass through.

Contraindications:

- Do not use in patients with known hypersensitivity to polyester.
- It should not be used in regions where it is not possible to obtain adequate adaptation of the band, compromising the proximal anatomy or cervical sealing.
- Contraindicated for deep subgingival procedures, in which the band cannot guarantee isolation, adaptation or stability.
- Do not use the damaged or deformed product.

Instructions for use:

- Wrap the tooth with the die;
- Adjustment to the desired size;
- For a better adaptation, use wedges in the interproximals of the tooth;
- Perform the restoration;
- Remove the wedges;
- Remove the matrix band.

Form of presentation

- Non-sterile product;
- One-time use
- Pack containing 25, 50 or 100 units.

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

- Non-sterile product.
- Single-use product.
- Do not subject the device to high temperatures.
- Before use, inspect that the Polyester Matrix Band is in perfect condition for use.

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REV. 01

- During operation, the professional must use all personal protective equipment, including goggles.
- Wear protective glasses on the patient.
- The Polyester Matrix Band is considered a contaminated instrument and proper handling precautions are required. Wear PPE, including gloves, masks, safety glasses, and a long-sleeved lab coat when handling it, as with other dental instruments.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Conservation and storage:

- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (Qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Polyester

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	It indicates a health product that has not been subjected to a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure
	Medical device	It indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directive 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.

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REV. 01

	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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