ENGLISH

Indication:

The carbide dental drill is a cutting device for use in dental treatments that, coupled to a high or low rotation micromotor, is used to act on a wide variety of tissues, such as: bone, enamel, dentin and restorative materials, such as hard metals and composites.

a) Special types of carbide drills

- Zekrya: Surgical use, odontosection, root separation, removal of root debris, the indication of its size varies according to the region to be reached.
- Endo-Z: endodontic use, preparation/smoothing of the pulp chamber during access surgery.
- Transmetal: High cutting power, it can remove restorative materials such as amalgam, in addition to the removal of metal prosthetic parts such as metal crowns, cast metal restoration, among others.
- Lindemann and Lindemann Cross Cut: Surgical use, indicated for cutting hard tissues/bone. The cross-cut ensures a more efficient cut.
- Surgical Drills: Used for surgical procedures, bone cuts, and dental cuts.
- Finishing Drills - 12 blades: Used in the finishing of restorations and rounding of cavo-superficial angles, the choice of the drill is preferably made by the professional performer according to the region of the oral cavity and dental element.
- Finishing Drills - 30 blades: Used in finer and more delicate finishes in restorations and rounding of the cavo-superficial angle, it can also be used for the removal of orthodontic resin remains. The choice of drills is preferably the responsibility of the professional performer according to the region in the oral cavity and dental element.

b) FG carbide drill bits

- Spherical: Used for cavity access, eventually for opening cavities in order to remove decayed tissue, delimiting additional retention and finishing line for prosthetic preparations.
- Pear: It is used to make additional retentions, rounding of angles, and eventually finishing lines and cingulate finishing.
- Inverted Conical: Cavity preparations where the walls must converge to the occlusal and create retentions.
- Cylindrical: Access to regions of greater difficulty, finishing of parallel surrounding walls, pulp wall planning.
- Conical: Preparation for indirect restorations, making grooves and grooves and eventually additional retentions.

c) CA carbide drill bits

These drills are used in a low rotation micromotor (contra-angle), the indications are the same as the FG, but with a lower rotation, ideal for the removal of decayed tissue near the pulp where the work must be meticulous to avoid its exposure. Despite the lower rotation, it is still necessary to use air/water cooling to avoid the heat generated by the friction of the drill.

Intended user:

- Duly qualified dental professionals.

Target population:

- The product is intended for use in dental patients undergoing restorative or orthodontic procedures (interproximal wear) as indicated by the dentist.

Fundamentals of its operation and action:

Carbide drills are mounted on low-speed or high-speed handpieces (dental turbines) that have high rotations and low torque, so they are small diameter grinding wheels. The activation of these instruments can be carried out by means of a compressed air turbine or an electric motor. Compressed air transmission generally provides lower available torque value and higher

DENTAL CARBIDE BUR
FRESA DE CARBURO DENTAL
BROCA CARBIDE ODONTOLÓGICA

REV.01

rotation. For electric motor drive, the available torque is higher; But the rotation is lower. Carbide Drills are produced in Tungsten Carbide and Stainless Steel, with measurement, tolerance, tensile strength and hardness.

Contraindications:

- To avoid fracture and deformation of the instrument, do not apply too much force when inserting the carbide drill bit into the handpiece (dental turbine).
- Respect the working speed (rpm Max), otherwise it may generate unwanted heat.
- Do not use the carbide drill bit as leverage and/or at incorrect angles.
- During use, do not apply excessive pressure to the carbide drill bit, as this may generate unwanted heat and/or cause carbide drill bit failure and risk to the patient.
- Prolonged immersion of carbide drill bits in disinfectant solutions can cause corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) for sterilization of the dental carbide drill. These chemical agents can cause corrosion of the products.
- Do not use sodium hypochlorite as a cleaning agent as it may cause corrosion in the product.
- Do not use the device on patients who are allergic to any tungsten carbide component.

Instructions for Use:

1. Insert the FG carbide bit into the high-speed turbine coupling using the knob or drill puller.
2. In the case of CA carbide drill bits, unlock the counter-angle fitting coupled to the low-speed micromotor, insert the drill bit rod and lock it properly.
3. Select the carbide drill bit according to the procedure performed and follow the speed according to the table below.
4. Check that the drill bit is correctly seated on the rotating tool and that the assembly works properly before use.
5. Bring the micromotor to the desired location.
6. Activate the micromotor.
7. Apply the active drill bit to the desired location intermittently.

HIGH ROTATION - FG CARBIDE

Diameter ISO	Speed Recommended	Maximum speed	Working pressure
(1/10 mm)	RPM (Rotação por minuto)	RPM (Rotação por Minuto))	Peso (g)
005-018	250.000 – 300.000	450.000	20 – 50
019-031	100.000 – 120.000	160.000	20 – 50
032-047	30.000 – 90.000	120.000	20 – 50
049-063	10.000 – 25.000	80.000	20 – 50
064-090	5.000 – 15.000	45.000	20 – 50

LOW ROTATION: CARBIDE CA

Diameter ISO	Speed Recommended	Speed Maximum	Pressure Work
(1/10 mm)	RPM (min ⁻¹)	RPM (min ⁻¹)	Peso (g)
005-014	70.000 – 90.000	250.000	30 – 50
015-023	60.000 – 80.000	150.000	30 – 50
024-040	40.000 – 50.000	80.000	50 – 100
041-100	20.000 – 40.000	50.000	50 – 100
110-220	7.000 – 12.000	20.000	50 – 100

Icon Automatic description

USE LIMITATIONS: Only 5 reuses are recommended.

Carbide drill bit cleaning can be performed in two ways: Manual or Ultrasonic. The professional who performs the cleaning and sterilization of the device must be using all protective PPE to avoid the risk of cross-contamination.

Manual Cleaning Procedure:

- Use an area specifically set aside for cleaning.

This review cancels and replaces the previous revision
3R Indústria e Comércio Ltda.
 Rua Ptolomeu, 290 – Vila Socorro
 São Paulo - SP – Brasil - CEP: 04762-040

MS: 80676920039

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DENTAL CARBIDE BUR
FRESA DE CARBURO DENTAL
BROCA CARBIDE ODONTOLÓGICA

REV.01

- Rinse the dental carbide drills under running water and keep them immersed in an enzymatic detergent solution bath, following the manufacturer's instructions.
- After immersing the device and keeping it submerged, brush it thoroughly with the neutral cleaning agent, directing the bur away from your body. Care should be taken to avoid spreading contaminants from the spray during brushing procedures. Do not use wire brushes.
- Special care should be taken to clean cracks and other hard-to-reach areas. Visually inspect to confirm debris removal, repeat cycle if necessary.
- Rinse the carbide drill bit under running water until it is visibly clean.
- Dry the appliance with clean compressed air or paper towels.

Ultrasonic Cleaning Procedure:

- Prepare an enzymatic detergent solution, following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place the carbide drill bit in the appropriate location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse the carbide drill bit under running water until it is visibly clean.
- Dry the appliance with clean compressed air or paper towels.

Sterilization requirements with carbide drills:

- This item applies to carbides drills purchased in non-sterile form and to the reuse of products purchased sterile.
- It is recommended that the carbide drill be sterilized by pressure saturated steam autoclave at 121 °C for 20 minutes, properly packed in its own autoclave packaging or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used.
- It is the professional's sole responsibility to keep the autoclave equipment up to date to ensure the correct sterilization of the product.

Warnings and Precautions:

- Before use, check that the Carbide Drill is in perfect condition for use, especially in relation to the integrity of the shank.
- Carbide drill bits are compatible with all handpieces that meet the requirements of applicable IEC 60601 series standards.
- Before use, check that the drill bit is securely seated in the low or high speed micromotor and that the coolant is working properly.
- Never force the carbide drill bit into the turbine. In case of mounting difficulty, check the turbine and the carbide dental drill.
- Respect the working speed (rpm Max), otherwise it may generate unwanted heat.
- Use proper cooling methods to ensure proper cooling of the parts used in the treatment.
- Use suction techniques that do not compromise the cooling of the parts involved in the treatment.
- Safety glasses should be worn to protect against ejected particles.
- Wear personal protective equipment, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the carbide drill bit, as you would with all other dental instruments. Wear protective glasses on the patient as well.
- In cases of returning the carbide drill bit, make sure it has been properly cleaned, sterilized, and packaged for return.
- BATCH NUMBER: It is important to keep the batch number, identified on the product label, to be informed of possible incidents with the manufacturer.

Discard:

- The disposal of the material must follow the relevant local legislation.

DENTAL CARBIDE BUR
FRESA DE CARBURO DENTAL
BROCA CARBIDE ODONTOLÓGICA

REV.01

Shelf Life:

- Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

- 5 uses (5 reuses), however, it is of paramount importance that before each use the check is carried out, that is, if the carbide drills are in perfect condition for use.

Conservation and Storage:

- This product is supplied non-sterile and must be cleaned and sterilized before first use and between uses, as described in these instructions for use.
- The product must be stored in a ventilated, dry place protected from light to maintain optimal storage and transport conditions, as well as physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed of the contents. The precautions that must be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packaging should not be left to the action of time.

Additional information:

- The sterilization processes described above have been validated and are suitable and effective for the use of carbide cutters. 3R is not responsible if the instructions for use are not followed correctly, generating unintended results for the efficacy and safety of the product/patient.
- If necessary, 3R provides instructions for use in printed format, at no additional cost, upon request by e-mail (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will return it to you within a maximum of 7 business days from the request.
- Always consult the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

- Tungsten carbide

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.

DENTAL CARBIDE BUR
FRESA DE CARBURO DENTAL
BROCA CARBIDE ODONTOLÓGICA

REV.01

	Maker	Indicates the manufacturer of the health product as defined in Directive 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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