BROQUEIRO DE ALUMÍNIO
ALUMINUM BUR BLOCK
FRESERO DE ALUMINIO

REV.01

ENGLISH

Indication:

The product is indicated to easily organize drills and diamond dental tips of various sizes and allows their sterilization.

Intended user:

Dental professionals duly qualified to manipulate, organize and select drills during clinical procedures.

Fundamentals of its operation and action:

Suitable for easily organizing drills and diamond dental tips of various sizes and allows sterilization.

Instructions for use:

- Make sure you have the correct size of the drill bit or diamond tip.
- Place the drill bits/bits in the corresponding holes once the size of the drill bits/bits has been determined.
- Close the Broqueiro.
- Proceed with the sterilization process.

Form of presentation

Model	Product	Specification
	FG	05 furos – FG 10 furos – FG 15 furos – FG 21 furos – FG 24 furos – FG 60 furos – FG
	RA/CA	05 furos – RA/CA 10 furos – RA/CA 15 furos – RA/CA 21 furos – RA/CA 24 furos – RA/CA 60 furos – RA/CA
	Broqueiro PM/HP	01 furos – HP/PM 05 furos – HP/PM 15 furos – HP/PM

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.

BROQUEIRO DE ALUMÍNIO
ALUMINUM BUR BLOCK
FRESERO DE ALUMINIO

REV.01

- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking, brush the device thoroughly away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by pressure saturated steam autoclave at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used. It is the professional's full responsibility to keep the autoclave equipment in ideal conditions of use to ensure the correct sterilization of the product.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Contraindications:

Not applicable

Warnings and Precautions:

- The cleaning of the Broqueiro can be done with a brush suitable for cleaning, water and neutral soap. Clean and dry material should follow autoclave sterilization procedures.
- The material must be packed in surgical-grade paper for sterilization.
- Never touch the Broqueiro to the wall of the autoclave, as they may be damaged.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Conservation and storage:

- Store in a clean, dry place.
- Transport in a clean, dry place at room temperature.

Discard:

The disposal of the material must follow the relevant local legislation.

BROQUEIRO DE ALUMÍNIO
ALUMINUM BUR BLOCK
FRESERO DE ALUMINIO

REV.01

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

The service life of the product may vary according to the conditions of use, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Aluminium

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	Date of manufacture	Indicates the date the health product was manufactured.
	Maker	Indicates the manufacturer of the health product as defined in Directive 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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