BROQUEIRO DE POLICARBONATO
POLYCARBONATE BUR BLOCK
FRESERO DE POLICARBONATO

REV.01

ENGLISH

Indication:

The product is indicated to easily organize dental burs of various sizes.

Intended user:

Dental professionals duly qualified to manipulate, organize and select burs during clinical procedures.

Fundamentals of its operation and action:

Suitable for easily organizing dental burs of various sizes.

Contraindications:

- Do not autoclave

Instructions for use:

- Make sure of the correct size of the bur.
- Place the burs in the corresponding holes once the size of the burs has been determined.
- Place the lid on the base, closing the bur block using the magnets.

Form of presentation

- Bur block CA – 35 holes
- Bur block FG – 35 holes
- Bur block PM – 7 and 21 holes

Manual Cleaning Procedure

To carry out manual cleaning, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking, brush the device thoroughly away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using clean compressed air or a paper towel.
- Repeat the cycle if necessary.

BROQUEIRO DE POLICARBONATO
POLYCARBONATE BUR BLOCK
FRESERO DE POLICARBONATO

REV.01

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using clean compressed air or a paper towel.

This device is not self-cleaning.

Warnings and Precautions:

- Do not use cold sterilization methods (chemical process) for the sterilization of the bur block. These chemical agents can cause corrosion of the products.
- Do not perform sterilization by autoclave.
- LOT Number: It is important to keep the lot number, identified on the product label, so that it can be reported to the manufacturer in the event of any occurrence.

Conservation and storage:

- Store in a clean, dry place.
- Transport in a clean, dry place at room temperature.

Disposal:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

The service life of the product may vary according to the conditions of use, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

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BROQUEIRO DE POLICARBONATO
POLYCARBONATE BUR BLOCK
FRESERO DE POLICARBONATO

REV.01

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the medical device should not be used.
	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	Date of manufacture	Indicates the date the medical device was manufactured.
	Manufacturer	Indicates the manufacturer of the medical device as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the medical device can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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