ENGLISH

Indication:

The EndoAccess kit diamond tip for endodontics is specifically designed for cavity access and preparation procedures in endodontic treatments. Its conical shape with rounded end allows the precise opening of the access cavity, in addition to favoring the uniform preparation of the pulp chamber walls. The model with an inactive tip offers greater safety to the dentist when working in the ceiling and floor regions of the pulp chamber, minimizing the risk of perforations and providing better control during wear.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who need endodontic treatment involving coronary opening, access to the pulp chamber and initial preparation of root canals, whether single-rooted or multi-rooted teeth.

Mechanism of action:

The diamond tips act by controlled abrasion, removing dental tissue precisely during access preparation. The conical shape and rounded end allow you to delimit the cavity, expand the visual field, and expose the pulp chamber safely.

Contraindications:

- To avoid fracture and distortion of the instrument, do not apply excessive force when inserting the diamond tip for endodontics into the high-rotation handpiece;
- Respect the working speed (RPM Max) defined on the product label, otherwise it may generate undesirable heat;
- Do not use the diamond tip for endodontics as levers and/or at incorrect angles;
- Do not apply excessive pressure to the diamond tip for endodontics, as this can cause undesirable heat and/or cause diamond tip failures leading to possible fracture.
- Cleaning agents with chlorine or chloride as an active ingredient are corrosive to stainless steel and should not be used.
- Prolonged immersion of diamond tips for endodontics in disinfectant solutions can result in corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) for sterilization of diamond tips for endodontics. These chemical agents can cause corrosion of the products.

Instructions for use:

- Insert the diamond tip for endodontics into a high-speed handpiece that meets the applicable requirements and standards for the product.
- Before activating the high-speed handpiece, check that the diamond tip is properly fitted and attached to the instrument.
- Activate the instrument and apply the active part of the diamond tip for endodontics to the surface that will be worn.
- As it is a diamond tip that will work at high speed, it is necessary that during use the professional is constantly using coolant.
- At the end of the procedure, sanitize the diamond tip and sterilize it through a saturated steam autoclave under pressure at 121°C for 20 minutes.

Select the diamond tip for endodontics that will be used and check the RPM according to the table below:

Diameter	Recommended Speed	Maximum Speed	Working pressure
----------	-------------------	---------------	------------------

ISO			
(1/10mm)	RPM (Rotation per Minute)	RPM (Rotation per Minute)	Weight (g)
008-018	250.000 – 300.000	450.000	20 – 50
021-031	100.000 – 120.000	160.000	20 – 50
033-047	30.000 – 90.000	120.000	20 – 50
049-063	10.000 – 25.000	80.000	20 – 50
065-090	5.000 – 15.000	45.000	20 - 50

Table: RPM (FG)

Form of presentation

- Non-sterile product;
- Kit with bur block + 5 diamond tips;

Cleaning procedures (diamond tip for endodontics):

- **Clean** and sterilize the diamond tip for endodontics before each use according to cleaning and sterilization item.
- It is recommended that cleaning be done immediately after the end of the procedure, to avoid drying out contaminants, which can make cleaning difficult.
- **After** cleaning, check the diamond tip for endodontics to ensure that any traces of contamination have been removed. Repeat cleaning if necessary.
- Cleaning the diamond tip for endodontics can be done in two ways: Manual or Ultrasonic.

Manual Cleaning Procedure:

If manual cleaning is the only option available, proceed as follows.

- Use an area specifically reserved for cleaning;
- Wash the diamond tip for endodontics under running water and keep them immersed in a bath of neutral pH cleaning solution (enzymatic detergent) following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent (enzymatic detergent). Care should be taken to avoid spreading contaminants by spraying during the brushing process. If necessary, it is possible to use specific brushes for cleaning dental drills.

- Care should be taken to clean cracks and other hard-to-reach areas. Visually inspect to confirm debris removal. Repeat the cycle if necessary.
- Rinse the endodontic diamond tip thoroughly under running water until they are visibly clean.
- Dry the device with clean compressed air or paper towels.

Ultrasonic Cleaning Procedure:

For ultrasonic cleaning, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent) following the agent manufacturer's instructions for correct concentration, exposure time, temperature, and water quality.
- Place the diamond tip for endodontics in the appropriate location for the ultrasonic equipment, following the manufacturer's instructions for use.
- Completely immerse the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse the endodontic diamond tip thoroughly under running water until they are visibly clean.
- Dry the device with clean compressed air or paper towels.

Sterilization requirements for diamond tip for endodontics:

- It is recommended to sterilize the diamond tip for endodontics in an autoclave with saturated steam under pressure at 121°C for 20 minutes, or to follow the sterilization parameters according to the relevant policies/regulations of the country where the device is being used.
- It is the professional's sole responsibility to keep the autoclave equipment up to date to ensure the correct sterilization of the product.

Warnings and Precautions:

- Before and after using the diamond tip for endodontics, proceed according to the cleaning and sterilization items.
- Clean and sterilize the diamond tips for endodontics before each use according to cleaning and sterilization items.
- For user and patient safety, the diamond tip for endodontics should only be used by qualified dental professionals.
- Before use, check that the diamond tip for endodontics is in perfect condition for use, especially in relation to the active area and integrity of the stem.
- Before use, make sure that the rod of the diamond tip for endodontics is fully fitted to the high-speed handpiece;
- Never force the diamond tip for endodontics into the high-speed handpiece. In case of difficulty in fitting, check the turbine and the diamond tip for endodontics.
- Respect the working speed (RPM Max) defined on the product label, otherwise you may generate undesirable heat.

Using:

- Move the diamond tip for endodontics intermittently to avoid localized heating and consequent breakage of the diamond tips, causing risks to the patient;
- Eye protection is of paramount importance to protect against ejected particles.
- Used diamond tips for endodontics should be considered contaminated instruments and appropriate handling precautions are required. Wear personal protective equipment (PPE), including gloves, masks, and safety glasses when handling, as with all other dental instruments.

Post-Use:

- Reuse of the product causes wear of the active part reducing its performance.
- In cases of return of the diamond tip for dental endodontics to 3R (manufacturer), make sure that they have been effectively cleaned, sterilized and properly packaged for return.

KIT ENDOACCESS
ENDOACCESS KIT
KIT ENDOACCESS

REV.01

- LOT Number: It is important to keep the lot number, identified on the product label, to be informed of any occurrences with the manufacturer.
- Material must be disposed of in accordance with local regulations.

Storage conditions:

- Storage should be dry and cleaned at room temperature.
- Keep the product stored in its original packaging until the moment of use.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

USE LIMITATIONS: The recommended shelf life of the product is 5 uses (5 reuses). However, it is of paramount importance that before each use, a check is carried out to ensure that the diamond tips are in perfect condition. Otherwise, they can overheat the dental tissues, cause a pulp response, and trigger a painful sensation in the patient, as well as not meeting the necessary wear requirement.

Additional Information

1. The reuse processes described above have been validated and are suitable and effective for the reuse of diamond tip for endodontics. 3R is not responsible if the instructions for use are not properly followed, resulting in undesirable results in product/patient efficacy and safety.

2. If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request via email (qualidade@3r.ind.br) or telephone (11-5525-0590 / 0800-773-5525), and we will get back to you no later than 7 business days from the patient's request.

Always check the website: www.microdont.com.br, the most current version of the usage instruction corresponding to the product.

Composition:

Active part: diamond

Stem: AISI 303 STEEL

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.

	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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