ENGLISH

Indication:

Prophylaxis brushes have a cup shape and as their name suggests they act in the prophylaxis and removal of bacterial plaque, in addition to having specific indications according to the thickness of their bristles/color, namely:

Pink – Extra Soft: Indicated for dental prophylaxis procedures, especially in patients with greater gingival sensitivity. Its extra-fine bristles promote gentle cleaning and minimize the risk of trauma to soft tissues.

Lilac – Soft: Indicated for dental prophylaxis with efficient removal of bacterial plaque. The thin bristles give it a soft characteristic, favoring patient comfort and reducing the risk of gingival trauma.

Yellow – Medium: Indicated for dental prophylaxis, with greater abrasive capacity compared to soft versions. Suitable for aid in the removal of extrinsic stains, as long as they are used with control and appropriate technique.

Blue – Semi-Rigid: Indications for dental prophylaxis and for clinical situations that require greater mechanical action, such as controlled microabrasion procedures of the enamel surface.

Intended user:

Duly qualified dental professionals.

Target population:

Patients who need dental prophylaxis procedures to remove bacterial plaque, maintain oral health and, when indicated, help remove extrinsic stains. It ranges from patients with greater gingival sensitivity, who benefit from extra-soft brushes, to those who require more intense mechanical action, such as in cases of microabrasion or the need for greater abrasiveness. The selection of the appropriate brush should consider the periodontal condition, the level of sensitivity and the specific needs of each patient.

Fundamentals of its operation and action:

The performance of brushes for prophylaxis is based on the mechanical friction between the bristles and the dental or restorative surfaces. The abrasive contact generated by the controlled rotational movement promotes the physical removal of dental biofilm, food residues and extrinsic stains, contributing to the finishing and clinical hygiene of the treated surfaces.

Contraindications:

- Regions with intense gingival inflammation, profuse bleeding or ulceration.
- Patients with known allergy to the components of this product.
- Presence of open carious lesions, extensive fractures or areas of severely weakened enamel.
- Subgingival areas, where use of the brush is not indicated due to the risk of tissue trauma.

Instructions for use:

The prophylaxis brush is a rotating instrument, used by dental professionals and mounted on a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor:

- Unlock the low-speed micromotor fitting;
- Insert the rod into the socket; Lock the rod on the low-rpm micromotor;
- Check that the rod is properly fitted to the low-speed micromotor and that the coolant is working properly;
- Place prophylaxis paste in the dappen dish and load the brush with the paste;
- Take the micromotor with the brush to the tooth surface;

ESCOVA PARA PROFILAXIA - COLORIDAS
PROPHYLAXIS BRUSHES - COLORED
CEPILLOS DE PROFILAXIS – COLOREADOS

REV.01

- Drive the micromotor at low speed without applying pressure so as not to generate excessive heat and move slowly evenly and humidify the tooth surface;
- If necessary, repeat the procedure.

Form of presentation

- Non-sterile product.
- Single use.
- Blue: Semi-rigid; Yellow: Medium; Lilac: Soft; Pink: Extra soft

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Ensure that the patient also wears protective goggles.
- Check that the rod is properly seated in the low-speed micromotor and that the coolant is working properly.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Disposal:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Additional information:

ESCOVA PARA PROFILAXIA - COLORIDAS
PROPHYLAXIS BRUSHES - COLORED
CEPILLOS DE PROFILAXIS – COLOREADOS

REV.01

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Metal rod
 Nylon bristles

Symbology:

Symbol	Title of Symbol	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	It indicates a health product that has not been subjected to a sterilization process.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Keep dry	Indicates that the product should be kept dry.
	See instructions for use	Indicates the user's need to consult the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured.
	Medical device	It indicates that it is a medical device.
	Maker	Indicates the Manufacturer of the product for health as defined in the European Community Directive 90/385/EEC, 93/42/EEC and 98/79/EC.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code for accessing the electronic version of the product's Instructions for Use (eIFU).

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