ENGLISH

Indication:

Fixing discs, saws, polishers, felt, strips and diamond tips for use in low speed motor or handpiece.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing dental procedures that require finishing, polishing, cutting or refinement of dental or restorative surfaces, performed with discs, saws, polishers, felts, strips or diamond tips coupled to a low-speed motor.

Fundamentals of its operation and action:

Fixation by means of screws, fitting or attached to discs, saws, sandpapers, polishers, felts or diamond tips, intended for use in low rotation motor or handpiece.

Contraindication:

- Do not use the device if it is damaged, deformed or with excessive wear, compromising the fixation or stability.
- Do not use on motors or handpieces incompatible with the device's attachment system.

Instructions for use:

- Make sure the low-rpm motor or handpiece is switched off before installation.
- Attach the chuck to the low-speed motor or handpiece, following the manufacturer's instructions.
- Attach the desired accessory (discs, saws, sandpaper, polishers, felts or diamond tips) according to its coupling system (screw, fitting or pressure), ensuring that it is firmly stabilized.
- Check the stability and alignment of the assembly before use.
- Start the low-speed motor and test the rotation without contact with the tooth to confirm attachment and absence of abnormal vibration.
- Start the procedure according to the indicated technique, maintaining light pressure, refrigeration when applicable, and continuous control of the operative field.

Form of presentation:

- Non-sterile product;
- Package containing one unit.

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning;
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the agent manufacturer's instructions.

- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent (enzymatic detergent). Care should be taken to avoid spreading contaminants by spraying during the brushing process. Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas.
- Visually inspect to confirm debris removal.
- Repeat the cycle if necessary.
- Rinse thoroughly under running water until visibly clean.
- Dry the device using a clean compressed air, or paper towel.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (5) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry the device using a clean compressed air, or paper towel

Sterilization Requirements

Sterilization in a saturated steam autoclave under pressure at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used.

Warnings and Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- During operation, the professional must use all personal protective equipment, including goggles.
- Wear protective glasses on the patient.
- Check that the rod is properly seated in the low-speed micromotor and that the coolant is working properly.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

The shelf life of the product may vary according to the conditions of use, frequency of sterilization, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

MANDRIL
MANDREL
MANDRIL

REV.01

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Metal: Brass or steel

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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