ENGLISH**Indication:**

Its function is to fill the porosities of the surfaces created by the previous conditioning with phosphoric acid and allow the union between the enamel and dentin with the restoration materials, such as composite resin and resinous cements in the cementation of fiberglass posts, prosthetic pieces (crowns, inlays/onlays, veneers).

Intended user:

Duly qualified dental professionals.

Target population:

Indicated for patients who require adhesive restorations, including enamel and dentin procedures involving composite resin, resin cements, or cementation of prosthetic posts and parts. It covers adults, adolescents and, when clinically appropriate, children, as long as the treatment requires the use of light-curing adhesive systems.

Fundamentals of its operation and action:

The Light-curing Adhesive is an adhesive for enamel and dentin, and acts as a bonding agent between the dental structure and restorative materials. The Light-curing Adhesive combines primer and adhesive in a single bottle, requiring acid conditioning prior to adhesive application. Among the main characteristics are the high bond strength, high chemical stability and resistance of the adhesive film generated after polymerization and the adequate penetration capacity of the monomers in the collagen fibers of demineralized wet dentin, contributing to greater longevity and safety.

Contraindication:

- The product is contraindicated for patients with allergies to any of the components of the formula.

Instructions for Use:

- Perform prophylaxis with oil-free prophylactic paste and isolation of the operative field prior to the adhesive steps.
- Do total conditioning with 37% phosphoric acid for 30 seconds on enamel and 15 seconds on dentin.
- Wash the area thoroughly and then dry it in enamel with an air jet and dentin with cotton so as not to dehydrate it.
- Apply a thin layer of Light-curing Adhesive, with the help of a microbrush in active friction movements. Apply a jet of air to volatilize the solvent and light-cure for 20 seconds.

Form of presentation:

- Pack of 4 ml.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As it is a chemical product, the adhesive should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

- Close the bottle immediately after use.
- Products such as eugenol or clove oil inhibit polymerization.
- Exclusive use of the dental surgeon.

MICROBOND – ADESIVO FOTOPOLIMERIZÁVEL
MICROBOND – LIGHT-CURED ADHESIVE
MICROBOND – ADHESIVO FOTOPOLIMERIZABLE

REV.01

- Do not use the product after the expiration date.
- Keep out of reach of children.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Disposal:

The disposal of the material must follow the relevant local legislation.

Expiration date:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

Single-use product, with a service life restricted to the moment of clinical application. After being removed from the original packaging and used, any remaining quantity should be discarded immediately.

Conservation and storage:

- Store the bottle in a dry, cool, dust-free place, avoiding heat, sunlight, and other sources of direct light.
- Flammable product.
- Do not store near materials that have eugenol or clove oil in their composition.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Monomeric mixture and Silanized Fumed Silica.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Consult the instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Temperature limit	Indicates the temperature limits to which the health product can be safely exposed.
	Care	It indicates the need for the user to refer to the instructions for use of important preventive uses, such as warnings and precautions that, for various reasons, cannot be displayed on the medical device itself.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.

MICROBOND – ADESIVO FOTOPOLIMERIZÁVEL
MICROBOND – LIGHT-CURED ADHESIVE
MICROBOND – ADHESIVO FOTOPOLIMERIZABLE

REV.01

Symbol	Symbol Title	Description
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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