ENGLISH

Indication:

MICROSIL is used as a bonding agent to increase adhesion between organic materials — such as adhesives that contain methacrylic monomers in their composition — and inorganic materials. It is widely used in restorative procedures involving materials such as composite resin, fiberglass pins, resinous cements, glass-phase ceramics, and ceromeres.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in patients undergoing dental restorative procedures involving indirect and, when applicable, direct restorative materials, in which it is necessary to promote or optimize adhesion between organic and inorganic surfaces.

Fundamentals of its operation and action:

The product consists of silane, a bifunctional molecule that acts as a binding agent between inorganic surfaces containing silica and resinous materials. After being applied to the treated surface, the silane reacts with hydroxyl groups present in the glass phase of the substrate, forming stable chemical bonds. Simultaneously, its methacrylic groups are able to copolymerize with the monomers present in adhesive systems and resin cements. This mechanism favors the formation of a stable adhesive interface, contributing to the increase of bond strength between materials.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use on surfaces or materials that do not have a glass phase or silica in their composition, such as polycrystalline ceramics (e.g., zirconia or alumina), in which silane does not have adequate adhesive effectiveness.
- Do not use the product for purposes other than those for which it was developed.

Instructions for Use:

Cementing of fiberglass pins:

- Wipe the surface of the fiberglass pin with 70% alcohol and air blast dry.
- Apply **SILANE** with the help of a brush or microbrush and wait for it to dry for 1 minute.
- Cement the pin according to the recommendations of the adhesive system and the selected resin cement.

Repairs of ceramic composites:

- Perform prophylaxis of the surface to be prepared with pumice stone and water. Then, dry and isolate the operative field.
- Remove friable edges and make a bevel of approximately 2 mm on the surface with diamond tips.
- Apply the porcelain conditioner to the beveled surface for the time recommended by the ceramic manufacturer.
- Wash thoroughly and dry completely.
- Apply **SILANE** with the help of a brush or microbrush to the treated surface and wait for it to dry for 1 minute. If necessary, apply a second layer and wait again for drying for 1 minute.
- Apply the selected adhesive system and perform the light-curing according to the manufacturer's instructions.
- Perform the repair with composite resin, adapting the material to the margins and avoiding overcontours.

Ceramic Repair/Ceramic Fragment Bonding:

- Perform prophylaxis of the surface to be prepared with pumice stone and water. Then dry and isolate the operative field.
- Remove friable edges and roughen the fractured surface with diamond tips.
- Apply the porcelain conditioner to the prepared surface for the time recommended by the ceramic manufacturer.

MICROSIL – SILANO**MICROSIL – SILANE****MICROSIL – SILANO**

REV.01

- Wash thoroughly and dry thoroughly.
- Apply **SILANE** with the help of a brush or microbrush to the entire treated surface and wait for it to dry for 1 minute. If necessary, apply a second layer and wait again for drying for 1 minute.
- Apply the selected conventional adhesive system and perform photocuring according to the manufacturer's instructions.
- Position the ceramic fragment, maintaining light pressure for adaptation, and perform photocuring according to the manufacturer's recommendations of the restorative material.

Cementing of indirect restorations with resin cements:

- Prepare the dental surface according to the technique recommended by the cementing system.
- Perform acid etching of the inner surface of the restoration using the porcelain conditioner for the time recommended by the manufacturer, according to the type of ceramic. Then wash and dry thoroughly.
- Apply **SILANE** with the help of a brush or microbrush and wait for it to dry for 1 minute. If necessary, apply a second layer and wait again for 1 minute.
- Cement the restoration using resin cement according to the product manufacturer's instructions.

Form of presentation

- Pack of 5 ml.

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

Not applicable. As a chemical, silane should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

- Immediately close the bottle after use to prevent evaporation. Keep in a dry place and away from light.
- Flammable product. Contains ethanol.
- Keep out of reach of children.
- Do not use the product with the expired expiration date.
- Non-sterile product.
- Product for exclusive use by dental professionals.
- The entire procedure should be performed, preferably, with absolute isolation of the operative field.
- Extensive repairs with composites should be performed in layers to reduce the effect of polymerization contraction.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

Single-use product, with a service life restricted to the moment of clinical application. After being removed from the original packaging and used, any amount remaining in the device should be discarded immediately.

Conservation and storage:

- Keep the product in its packaging always closed.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- Store in a ventilated and dry place, away from direct sunlight and at room temperature.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Dynasylan Memo (SILANO)	1,5 %
Glacial Acetic Acid	2,5 %
Ethyl Alcohol Absolute	94,0 %
Purified water	2,0 %

Symbology:

Symbol	Symbol Title	Description
	Expiry date	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Caution	It indicates the need for the user to refer to the instructions for use of important preventive uses, such as warnings and precautions that, for various reasons, cannot be displayed on the medical device itself.
	Flammable product	Indicates that the product contains flammable substances and should be kept away of ignition sources.
	Recyclable packaging	Indicates that the packaging is recyclable.
	Keep away from sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or subject to prescription, in accordance with applicable legislation.

MICROSIL – SILANO
MICROSIL – SILANE
MICROSIL – SILANO

REV.01

Symbol	Symbol Title	Description
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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