MIR - IPR SYSTEM
MIR - IPR SYSTEM
MIR - SISTEMA IPR

REV.01

ENGLISH

Indication:

MIR is mainly indicated for orthodontic strippings (a technique in which you want to create interproximal space through enamel wear) and can also be used for interproximal finishing and for removing excess cement and restorative materials. It has a sequence of six arcs, one saw and five abrasive strips ranging from extra-fine to extra-coarse for the MIR 1.0 and five arcs, one saw and four abrasive strips ranging from extra-fine to coarse for the MIR 2.0.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative or orthodontic procedures (interproximal wear) as indicated by the dentist.

Fundamentals of its operation and action:

Abrasion

Contraindications:

- Do not use in patients with known hypersensitivity to metals, especially stainless steel;
- Do not use in areas with severe inflammation or active bleeding;
- Do not use the product damaged, deformed or with loss of abrasive coating.

Instructions for use:

The MIR was developed to assist in the interproximal opening of the dental elements for orthodontic stripping. Each saw bow or abrasive strip has a numbering (in millimeters) corresponding to the size of the opening, to obtain this opening it must be used following the ascending order of arcs:

MIR 1.0:

- White – Serrated – 0.05mm
- Yellow – Extra Fine – 0.10mm
- Red – Fine – 0.20mm
- Blue – Regular – 0.30mm
- Green – Coarse – 0.40mm
- Purple – Extra Coarse – 0.50mm
- The arcs of the MIR 1.0 can be used freehand, with the aid of the handle or in a contra-angle.

MIR 2.0:

- White – Serrated – 0.05mm
- Yellow – Extra Fine – 0.05mm single sided and 0.10mm double sided
- Red – Fine – 0.10mm single sided and 0.20mm double sided
- Blue – Regular – 0.15mm single face and 0.30mm double face
- Green – Coarse – 0.20mm single sided and 0.40mm double sided
- The arcs of MIR 2.0 are exclusively manual.

MIR PRIME / MIR 1.0 – With Cable

Sequence of Use:


MIR - IPR SYSTEM

MIR - IPR SYSTEM


MIR - SISTEMA IPR

REV.01

MIR 2.0 – MANUAL

Sequence of Use:

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

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MIR - IPR SYSTEM
MIR - IPR SYSTEM
MIR - SISTEMA IPR

REV.01

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

- Before use, inspect if the product is in perfect condition of use;
- Non-sterile product.
- During operation, the professional must use all personal protective equipment.
- Wear protective glasses on the patient.
- The Product is considered a contaminated instrument and proper handling precautions are required. Wear PPE, including gloves, masks, safety glasses, and a long-sleeved lab coat when handling, as is done with all other dental instruments.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- Keep the packaging in a dry place, away from light and away from the action of bad weather. Keep the products in their original packaging until the moment of use.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- 3R is not responsible if the instructions for use are not followed correctly, generating undesired results in the efficacy and safety of the product/patient.
- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

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Polycarbonate

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product as defined in Regulation (EU) 2017/745 (MDR).
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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