PASTA DIAMANTADA
DIAMOND PASTE
PASTA DE DIAMANTE

REV.01

ENGLISH

Indication:

Diamond paste is made up of natural diamond powder, added to a pasty and water-soluble vehicle. Indicated for polishing restorations made in composite resins and ceramics. Offers shine and smoothness, providing greater naturalness in aesthetic restorations. It can be used by dental surgeons (intra-oral) and dental prosthesis technicians (extra-oral).

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing direct or indirect restorative procedures in composite resin or ceramics, who require finishing and polishing to optimize the aesthetics, comfort and durability of rehabilitations.

Fundamentals of its operation and action:

The Diamond Paste has the right consistency to be placed directly on the restored teeth without running into the oral cavity. The product can be applied with Robinson brush, rubber bowls for prophylaxis and felt discs. It is abrasive enough to efficiently remove all types of residues from the tooth surface and polish the restoration without causing undue abrasion of the tooth surface.

Contraindication:

The product does not present contraindications when used for the purpose for which it is intended, following its instructions for use.

Instructions for use:

- After initial finishing with cups and/or abrasive brushes specific for finishing and dental sculpture, apply the desired amount of diamond paste on the surface to be polished. Impregnate with the diamond paste, a felt disc adapted to a low-speed mandrel or a non-abrasive rubber bowl.
- Always start polishing at low speed, with circular and intermittent movements, in order to minimize heat generation.
- The diamond paste is presented in three different granulations, syringes with 3 or 6 grams, available for medium, fine and extra-fine coating particle size: 6, 3 and 1 μ .
- It is recommended to start polishing with the highest grain paste (6 μ), carefully clean the treated area, and finish with the lowest grain paste.
- Wash and dry surfaces.
- Reapply the paste whenever you deem necessary in order to produce smoother surfaces.

Form of presentation

- Product in the form of a paste;
- Syringes with 3 g or 6 g;
- Available in 1 μ m, 3 μ m and 6 μ m granulations;

Cleaning Procedure:

Not applicable. It is a chemical product and, therefore, should not undergo any type of cleaning procedure.

Sterilization Procedure:

PASTA DIAMANTADA
DIAMOND PASTE
PASTA DE DIAMANTE

REV.01

Not applicable. As it is a chemical product, diamond paste should not be sterilized or subjected to any decontamination method.

Warnings and Precautions:

- For dental use only.
- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product after the expiration date or if the packaging is violated.
- Keep out of reach of children.
- Do not ingest. If ingested, it may cause nausea.
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- Diamond paste should be stored in a clean, moisture-free environment.
- Do not store at temperatures above 49°C (120°F) or below 5°C.
- Keep the product in its packaging always closed.
- Avoid contact with moisture
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf life:

The product has a shelf life of 2 years, and must be used within this period to ensure its effectiveness and safety.

Service life:

The product's service life corresponds to its shelf life of 2 years, provided the packaging remains closed and the storage conditions described above are observed. Once the syringe has been opened, the paste must be used within the validity period; any amount dispensed for a procedure and not used must be discarded immediately and must not be returned to the packaging.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

PASTA DIAMANTADA
DIAMOND PASTE
PASTA DE DIAMANTE

REV.01

Composition:

Polyethylene glycol, glycerin, preservatives Methylparaben and propylparaben, and micronized diamond powder.

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Medical device	Indicates that it is a medical device.
	CE Mark	It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Maker	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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