POLIDORES PARA AMÁLGAMA
AMALGAM POLISHERS
PULIDORES DE AMALGAMA

REV.01

ENGLISH

Indication:

Amalgam polishers are suitable for finishing and polishing amalgam restorations. Its use provides smooth, uniform surfaces with adequate shine, contributing to better clinical performance and longevity of the restoration. The instrument must be coupled to a low-speed, compressed air-driven micromotor or electric motor.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the finishing and polishing of amalgam restorations, according to the clinical indication of the dental surgeon

Fundamentals of its operation and action:

Amalgam polishers work by means of controlled abrasive action. During rotation, the abrasive particles from the surface of the polisher come into contact with the restoration, promoting the removal of irregularities and the progressive smoothing of the amalgam. This process reduces surface roughness, improves gloss, and minimizes biofilm buildup, contributing to the durability and clinical performance of the restoration.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made of materials other than Amalgam.
- Do not use the product for purposes other than those for which it was developed.

Instructions for use:

Start polishing the metal with the most abrasive piece and work your way up to the least abrasive piece (in the following order: Coarse, medium, fine)

- Unlock the micromotor fitting
- Insert the rod into the micromotor
- Lock the rod on the micromotor
- Check that the rod is properly seated in the micromotor and that the coolant is working properly
- Apply polishing paste to the polishing piece.
- Bring the micromotor with the polishing part to the tooth surface
- Drive the micromotor at low speed without applying pressure so as not to generate excessive heat and move slowly evenly
- If necessary, repeat the procedure

Form of presentation:

- Non-sterile product;
- Models: Cup and torpedo;
- Kit with 6 polishers; and 8 polishers;
- 4 granulations: Coarse (Brown), medium (Green), Fine (Blue) and Extra fine (Beige).

POLIDORES PARA AMÁLGAMA
AMALGAM POLISHERS
PULIDORES DE AMALGAMA

REV.01

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by autoclave of saturated steam under pressure is recommended, at 121°C for 10 minutes, or sterilization parameters should be followed according to the relevant policies/regulations of the country where the device is being used.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Warnings and Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- Before use, inspect that the polishes are in perfect condition for use;
- Non-sterile product.
- Sterilize the product before use.
- Make sure that the rubber bowl/cone is properly seated in the low-speed motor.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Wear protective goggles on the patient.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

POLIDORES PARA AMÁLGAMA
AMALGAM POLISHERS
PULIDORES DE AMALGAMA

REV.01

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

5 uses (5 reuses), however, it is of paramount importance that before each use a check is carried out, that is, if the amalgam polishers are in perfect condition for use

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Silicone

Abrasive: Aluminum Oxide

Metal rod (Brass Ø 1/8 + nickel plating)

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Maker	Indicates the manufacturer of the health product as defined in Directive 90/385EED, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).

POLIDORES PARA AMÁLGAMA
AMALGAM POLISHERS
PULIDORES DE AMALGAMA

REV.01

	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).
---	--	--

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