ENGLISH

Indication:

The parts for resin polishing (Disc, torpedo, cup, torpedo) are rotating instruments, which provide greater smoothness and shine of resins, used by dental professionals, installed in a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor.

3R is not responsible for improper use of the product, especially when used for purposes other than those for which it was specifically developed.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the finishing and polishing of composite resin restorations, as indicated by the dentist.

Fundamentals of its operation and action:

Abrasion. Dental polishers have their action in polishing and finishing, they have varied formats for application on the composite resin on all faces.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than the composite resin.
- Do not use the product for purposes other than those for which it was developed.
- Do not use continuously; the product must be used with intermittent movements and with cooling, as excessive heat may cause damage to the dental tissue.

Instructions for use:

The parts for resin polishing (Disc, torpedo, cup, torpedo) are rotating instruments, which provide greater smoothness and shine of resins, used by dental professionals, installed in a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor.

- Carefully open the resin kit making sure all parts are inside the package.
- Immediately afterwards, attach the polisher shank to the low-speed handpiece, making sure it is well locked.
- Start the resin polishing procedure by applying a small amount of polishing paste to the surface of the resin.
- Soon after, with the tip properly fitted into the low rotation motor, using the thickest polishing tip on the surface of the resin, providing smoothness and shine.
- It is recommended to change from the thickest tip to the thinnest tip, (Coarse, Medium, Fine) without applying pressure, only lightly on the tooth surface, and without applying excessive lateral pressure.
- It is essential to use a polishing paste, for better performance of the tips.
- If necessary, repeat the procedure more than once.
- The recommended rotation range per minute is from 5,000 to 8,000 RPM, and you should not exceed 12,000 RPM under any circumstances.

POLIDORES PARA RESINA
RESIN POLISHERS
PULIDORES DE RESINA

REV. 02

Note: Some polishers have a CA or FG shank. The CA shank has a low-speed fitting, and the FG shank is a different shank, to be fitted at high speed or used with an adapter for low speed.

Form of presentation:

- Non-sterile product;
- Models: Torpedo, torpedo, cup and disc;
- Kit with 6 polishers; 8 polishers or 12 polishers;
- 3 granulations: Coarse (Green), Medium (Yellow) and Fine (White).

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by pressure saturated steam autoclave at 121 °C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used. It is the professional's full responsibility to keep the autoclave equipment in ideal conditions of use to ensure the correct sterilization of the product.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Warnings and Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- Before use, inspect that dental polishes are in perfect condition for use;

POLIDORES PARA RESINA
RESIN POLISHERS
PULIDORES DE RESINA

REV. 02

- Non-sterile product.
- Sterilize the product before use.
- Make sure the dental polishes are properly seated in the low-speed motor.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Wear protective goggles on the patient.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

Disposal:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

5 uses (5 reuses), however, it is of paramount importance that before each use a check is carried out, that is, if the dental polishers are in perfect condition for use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (Qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Silicone

Abrasive: Aluminum oxide ALL 320/600/1000

Metal rod (Brass Ø 1/8 + nickel plating)

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.

POLIDORES PARA RESINA
RESIN POLISHERS
PULIDORES DE RESINA

REV. 02

	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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