POLISOFT - POLIDOR DE FELTRO
POLISOFT - FELT POLISHER
POLISOFT - PULIDORA DE FIELTRO

REV.01

ENGLISH

Indication:

The impregnated felt polisher is a rotating instrument, which provides greater shine on the restored composite resin surface. As its name suggests, it has abrasive impregnated in its active part, which when in contact with water forms a paste on the restored tooth surface. It should be used in the final stage of finishing and polishing, installed in a contra-angle with a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor.

Fundamentals of its operation and action:

Abrasion. The impregnated felt polisher has its action in the polishing and finishing stage.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the finishing and polishing of composite resin restorations, as indicated by the dentist.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than composite resin.
- Do not use the product for purposes other than those for which it was developed.

Instructions for use:

The impregnated felt polisher is a rotating instrument, which provides greater smoothness and shine of composite resins. Used by dental professionals, installed in a contra-angle with a low-speed micromotor. The activation of this instrument can be done through compressed air or electric motor.

- Carefully open the package making sure that the piece is in perfect condition for use.
- Then install the felt disc on the mandrel and attach the mandrel rod to the contra-angle + low-speed micromotor assembly, making sure it is well locked.
- Start the procedure by humidifying the felt; to do this, apply a few drops of water on it and also on the tooth surface, this step is essential for the product to perform its function properly, as the abrasive will be activated and will form the paste for polishing.
- It is recommended that cooling be carried out between intermittent uses and that no pressure be applied to the tooth surface.
- If necessary, repeat the procedure more than once.

Form of presentation

- Non-sterile product;
- Single-use product;

The reprocessing of the felt disc is prohibited, that is, it must be discarded after use. The cleaning and sterilization procedures described below are only for the mandrel that comes with the product.

Manual Mandrel Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.

POLISOFT - POLIDOR DE FELTRO
POLISOFT - FELT POLISHER
POLISOFT - PULIDORA DE FIELTRO

REV.01

- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using clean compressed air or a paper towel.
- Repeat the cycle if necessary.

Ultrasonic Mandrel Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using clean compressed air or a paper towel.

Sterilization Requirements

- Sterilization by pressure saturated steam autoclave at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used. It is the professional's full responsibility to keep the autoclave equipment in ideal conditions of use to ensure the correct sterilization of the product.
- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Warnings and Precautions:

- Only legally qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.
- Before use, inspect whether the impregnated felt polisher is in perfect condition for use;
- Non-sterile product.
- Sterilize the product before use.
- Make sure that the felt polisher is properly seated in the low-speed motor.
- The product should be used between 4,000 and 10,000 RPM.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Also use protective goggles on the patient.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.
- Do not use in patients who have allergic reactions to nickel or other components present in the product.
- After use, the felt disc must be discarded, only the mandrel can be reprocessed.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

POLISOFT - POLIDOR DE FELTRO
POLISOFT - FELT POLISHER
POLISOFT - PULIDORA DE FIELTRO

REV.01

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

- Not applicable. The product is for single use only and must be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Active part – felt with aluminium oxide-based abrasive
 Stem – Nickel-plated brass

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the medical device should not be used.
	Non-sterile	Indicates a medical device that has not undergone a sterilization process.
	Keep dry	Indicates a medical device that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the medical device was manufactured.
	Keep out of sunlight	Indicates a medical device that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the medical device.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the medical device can be identified.

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