ENGLISH

Indication:

Robinson brushes are indicated for dental prophylaxis procedures, aiming at the removal of biofilm, extrinsic stains and debris from tooth surfaces. In the mini format, its reduced size and the concentration of bristles in a single active tip allow localized, precise application in hard-to-reach areas, such as cervical and interproximal regions and around restorations, prostheses and orthodontic accessories.

Intended user:

Duly qualified dental professionals.

Target population:

Robinson brushes are intended for patients of any age group who require professional prophylaxis, including biofilm removal, extrinsic stains, or cleaning hard-to-reach tooth surfaces.

Fundamentals of its operation and action:

Robinson brushes act by rotational mechanical action, promoting the removal of biofilm and surface pigments through controlled friction between their bristles and the tooth surface. This abrasive movement, associated with the use of prophylactic pastes, cleans and smooths the tooth surfaces, favoring oral hygiene.

Contraindications:

- Regions with severe gingival inflammation, profuse bleeding or ulcerations.
- Patients with a known allergy to the components of this product.
- Presence of open carious lesions, extensive fractures or areas with very fragile enamel.
- Subgingival areas, where the use of the brush is not indicated due to the risk of tissue trauma.

Instructions for use:

The Robinson Microdont Brush is a rotating instrument, used by dental professionals installed on a low rotation micromotor. The activation of this instrument can be done through compressed air or electric motor:

- Unlock the low-speed micromotor fitting;
- Insert the rod into the socket; Lock the rod on the low-rpm micromotor;
- Check that the rod is properly fitted to the low-speed micromotor and that the coolant is working properly;
- Apply prophylactic paste to the brush;
- Take the micromotor with the brush to the tooth surface;
- Drive the micromotor at low speed without applying pressure so as not to generate excessive heat and move slowly evenly;
- If necessary, repeat the procedure.

Form of presentation

- Non-sterile product.
- Single use.
- The traditional brushes are available in two active-tip shapes: conical and flat.
- The mini brushes are available in two versions, differentiated by the color and stiffness of the bristles:
 - White – soft bristles;
 - Black – semi-rigid bristles.

Cleaning Procedure:

Not applicable. The product is for single use, so it should not be cleaned or reused.

Sterilization procedure

ROBINSON - ESCOVA PARA PROFILAXIA
ROBINSON - PROPHYLAXIS BRUSHES
ROBINSON - CEPILLOS DE PROFILAXIS

REV. 01

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

Warnings and Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Wear protective goggles on the patient.
- Check that the rod is properly seated in the low-speed micromotor and that the coolant is working properly.
- Non-sterile product.
- Single-use product.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

Not applicable. The product is for single use only and should be discarded after use.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Additional information:

- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Metal rod
Nylon bristles

Symbology:

Symbol	Title of Symbol	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	It indicates a health product that has not been subjected to a sterilization process.

ROBINSON - ESCOVA PARA PROFILAXIA
ROBINSON - PROPHYLAXIS BRUSHES
ROBINSON - CEPILLOS DE PROFILAXIS

REV. 01

Symbol	Title of Symbol	Description
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Keep dry	Indicates that the product should be kept dry.
	See instructions for use	Indicates the user's need to consult the instructions for use.
	Date of manufacture	Indicates the date when the health product was manufactured.
	Medical device	It indicates that it is a medical device.
	Maker	Indicates the Manufacturer of the product for health as defined in the European Community Directive , 93/42/EEC and 98/79/EC.
	Lot Code	Indicates the manufacturer's batch code so that the lot or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code for accessing the electronic version of the product's Instructions for Use (eIFU).

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