ENGLISH**Indication:**

Sanding discs are suitable for finishing and polishing composite restorations and are designed to provide a smooth surface. For this, they should be used initially with the coarser granulation and then move on to the smaller granulations. They act as auxiliaries in the finishing and polishing stage, which must be complemented with the use of other polishers.

The discs are available in four grain sizes: coarse (black), medium (brown), fine (orange) and extra-fine (white). In addition, they are available in two sizes: Ø 14 mm and Ø 12 mm.

3R is not responsible for improper use of the product, especially when used for purposes other than those for which it was specifically developed.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in patients who have undergone restorative procedures, with the aim of providing the proper finishing and polishing of the restoration.

Fundamentals of its operation and action:

The polishing discs come with the pop-on mandrel, which provides versatility in use, as it is possible to easily remove the disc from the mandrel and invert it. This allows the disc to be applied to different regions of the tooth, since only one face of the disc is active, i.e. the abrasive is only present on one side.

Applications:

Finishing and polishing of restorations made in composite resin.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than the composite resin.
- Do not use the product for purposes other than those for which it was developed.

Instructions for Use:

- Fit the disc into the mandrel and make sure it is securely fastened.
- Fit the metal rod of the mandrel to the low-speed handpiece (contra-angle).
- Make sure that the metal rod and disc are properly seated in the low-speed motor.
- Adjust the speed to 8,000 rpm, do not exceed 10,000 rpm.
- During use, the professional must perform intermittent movements to avoid overheating and excessive wear.
- It is important to moisten the surface during use and wash thoroughly at grain changes in order to prevent granules from the previously used disc from scratching the surface when the smaller grain disc is applied.
- Discard it at the end of use.

The sequence of the procedure for using the discs should follow a decreasing order of granulation. Thus, the professional should start with the disc with the highest abrasiveness, the black one, and then progressively change to discs with less abrasiveness, in the order: black, brown, orange and, finally, white.

Form of presentation:

SUPERFLEX - POLIDONT

- Non-sterile product.
- Single use.
- Pack with 32 discs or refill with 8 discs.
- Black: Coarse / Brown: Medium / Orange: Fine / White: Extra-fine.

Cleaning Procedure:

- The discs are single-use; therefore, they should not be cleaned or reused.
- Only the mandrel can be cleaned and reused.

Sterilization Procedure:

- The discs are supplied non-sterile and should not be sterilized. They should be disposed of after use, according to local regulations for dental waste.
- Only the mandrel can be sterilized and reused.

Warnings and Precautions:

- To avoid fractures and distortion of the instrument, do not apply excessive force when inserting the disc into the low-speed handpiece.
- Before use, inspect if the disc is in perfect condition for use, especially in relation to the area with granulation and the integrity of the rod;
- Before use, make sure the disc rod is perfectly fitted to the contra-angle;
- Never force the disc into the contra-angle. In case of difficulty in fitting, check the contra-angle and the rod of the disc.
- It is recommended that the use of the product is always accompanied by absolute isolation of the operative field.
- Do not use the disc as a lever and/or at incorrect angles;
- Contact with soft tissues (gums, mucosa) can cause local lesions.
- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and should not be used.
- Prolonged immersion of the Disc in disinfectant solutions may result in corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) for disc sterilization. These chemical agents can cause corrosion of the products.
- Protective eyewear should be worn to protect against particles ejected on the patient and the professional.
- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- In case of returning the product, make sure that it has been effectively cleaned and properly packaged for return.
- **BATCH NUMBER:** It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

- The disposal of the material must follow the relevant local legislation.

SUPERFLEX - POLIDONT

REV. 01

Shelf life:

- Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

- Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (Qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Mandrel – AISI 304 stainless steel

Active Part - PET, polyester adhesive and alumina abrasive

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates a health product intended for a single use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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