ENGLISH

Indication:

Sanding discs are suitable for finishing and polishing composite restorations and are designed to provide a smooth surface. To do this, they should be used initially with the coarsest granulation and then with the finer granulations. They act as auxiliaries in the finishing and polishing stage, which must be complemented with the use of other polishers. The discs are available in four granulations identified by the color of the center: coarse (burgundy), medium (brown), fine (orange) and extra-fine (yellow). In addition, they are available in two sizes: Ø12.5mm and Ø9.5mm. The non-metallic center prevents possible scratches on the surface of the resin if this region accidentally touches the restored surface.

3R is not responsible for improper use of the product, especially when used for purposes other than those for which it was specifically developed.

Intended user:

- Duly qualified dental professionals.

Target population:

- The product is intended for use in patients who have undergone restorative procedures, with the aim of providing the proper finishing and polishing of the restoration.

Fundamentals of its operation and action:

The polishing discs come with the chuck that provides versatility in use, as it is possible to remove the disc from the chuck and invert it. This allows the disc to be applied to different regions of the tooth, since only one face of the disc is active, i.e. the abrasive is only present on one side.

Applications:

- Finishing and polishing of restorations made in composite resin.

Contraindications:

- The product should not be used in patients with a history of hypersensitivity to any components present in its composition.
- Do not use in restorations made with materials other than the composite resin.
- Do not use the product for purposes other than those for which it was developed.

Instructions for Use:

- Fit the disc into the chuck and make sure it is securely fastened.
- Fit the metal chuck shank to the low-rotation handpiece (counter angle).
- Make sure that the metal rod and disc are properly seated in the low-rpm motor.
- Adjust the speed to 8,000rpm, do not exceed 10,000rpm.
- During use, the professional must perform intermittent movements to avoid overheating and excessive wear.
- It is important to moisten the surface during use and wash thoroughly at grain changes in order to prevent granules from the previously used disc from scratching the surface when the lower grain disc is applied.
- Discard it when it is finished.

SUPERSOFT - POLIDONT
SUPERSOFT - POLIDONT
SUPERSOFT - POLIDONT

REV.

The sequence of the procedure for using the discs should follow a decreasing order of granulation. Thus, the professional should start with the disc with the greatest abrasiveness, the burgundy one, and then progressively change to discs with less abrasiveness, in the order: burgundy, brown, orange and, finally, yellow.

Form of presentation:

- Non-sterile product.
- Single use.
- Pack with 32 discs or refill with 8 discs.
- Burgundy: Thick / Brown: Medium / Orange: Fine / Yellow: Extra-fine.

Cleaning Procedure:

- The discs are single-use; therefore, they should not be cleaned or reused.
- Only the chuck can be cleaned and reused.

Sterilization Procedure:

- The discs are supplied non-sterile and should not be sterilized. They should be disposed of after use, according to local regulations for dental waste.
- Only the chuck can be sterilized and reused.

Warnings and Precautions:

- These discs are for single use only and are discarded after use. Only the mandrel can be sterilized and reused.
- To avoid fractures and distortion of the instrument, do not apply excessive force when inserting the disc into the low-speed handpiece.
- Before use, inspect if the disc is in perfect condition for use, especially in relation to the area with granulation and the integrity of the rod;
- Before use, make sure the disc rod is perfectly fitted to the contra-angle;
- Never force the disc into the turbine. In case of difficulty in fitting, check the turbine and the rod of the Disc.
- It is recommended that the use of the product is always accompanied by absolute isolation of the operative field.
- Do not use the disc as levers and/or at incorrect angles;
- Contact with soft tissues (gums, mucosa) can cause local lesions.
- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and should not be used.
- Prolonged immersion of the Disc in disinfectant solutions may result in corrosion and should therefore be avoided.
- Do not use cold sterilization methods (chemical process) for disc sterilization. These chemical agents can cause corrosion of the products.
- Protective eyewear should be worn to protect against particles ejected on the patient and the professional.
- Wear PPE, including gloves, masks, safety goggles, and long-sleeved lab coats, as is done with all other dental instruments.
- Wear protective glasses on the patient as well.
- In case of returning the product, make sure that it has been effectively cleaned and properly packaged for return.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

SUPERSOFT - POLIDONT
SUPERSOFT - POLIDONT
SUPERSOFT - POLIDONT

REV.

- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

- Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

- Not applicable. The product is for single use only and should be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (Qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Mandrel – AISI 303 stainless steel

Active part – PET and aluminum oxide abrasive

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be

SUPERSOFT - POLIDONT
 SUPERSOFT - POLIDONT
 SUPERSOFT - POLIDONT

REV.

Symbol	Symbol Title	Description
		identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code - Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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