TAÇA E CONE DE BORRACHA PARA ODONTOLOGIA
RUBBER CUP AND CONE FOR DENTISTRY
COPA Y CONO DE GOMA PARA ODONTOLOGÍA

REV. 01

ENGLISH

Indication:

The rubber cup is indicated for the mechanical removal of biofilm, extrinsic stains and softened residues from the tooth surface during dental prophylaxis procedures.

Intended user:

Duly qualified dental professionals.

Target population:

The rubber cup for prophylaxis is intended for patients of any age group who need professional removal of biofilm, extrinsic stains or dental polishing.

Fundamentals of its operation and action:

The rubber cup acts by rotational movement at low rotation, whose flexibility allows adaptation to dental surfaces, promoting mechanical removal of biofilm and extrinsic stains.

Contraindications:

- Regions with severe gingival inflammation, profuse bleeding or ulcerations.
- Patients with a known allergy to the components of this product.
- Presence of open carious lesions, extensive fractures or areas with very fragile enamel.
- Subgingival areas, where the use of the cup is not indicated due to the risk of tissue trauma.

Instructions for use:

- Prophylactic cups are rotating instruments installed on a low-speed micromotor used by dental professionals.
- The drive can be done by compressed air or an electric motor.
- Unlock the micromotor socket.
- Insert the rod into the socket.
- Lock the rod into the micromotor.
- Check that the rod is properly fitted and that the coolant is working properly.
- Apply prophylactic paste to the cup.
- Bring the assembly to the tooth surface.
- Activate the micromotor at low speed, without applying excessive pressure, moving slowly and evenly to avoid heating.
- Repeat the procedure if necessary.

Form of presentation

- Non-sterile product;

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush away from the body using the neutral cleaning agent (enzymatic detergent).
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.

TAÇA E CONE DE BORRACHA PARA ODONTOLOGIA
RUBBER CUP AND CONE FOR DENTISTRY
COPA Y CONO DE GOMA PARA ODONTOLOGÍA

REV. 01

- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by saturated steam autoclave under pressure is recommended, at 121 °C for 20 minutes, or follow the parameters required by local standards and regulations.
- Remove the product from the original packaging, pack it in a sterilization packaging and sterilize in an autoclave according to the instructions of the equipment.

Warnings and Precautions:

- Non-sterile product.
- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Wear protective goggles on the patient.
- Check that the rod is properly seated in the low-speed micromotor and that the coolant is working properly.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

Not applicable. The product is for single use only and should be discarded after use.

Additional information:

TAÇA E CONE DE BORRACHA PARA ODONTOLOGIA
RUBBER CUP AND CONE FOR DENTISTRY
COPA Y CONO DE GOMA PARA ODONTOLOGÍA

REV. 01

- If necessary, 3R provides instructions for use in printed format, at no additional cost, upon request by email (qualidade@3r.ind.br) or telephone (+55 11 5525-0590 / 0800 773 5525). The service will be carried out within 7 business days after the request.
- Always check the company's website for the most up-to-date version of the instructions for use.

Composition:

Natural silicone

Metal rod

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
Rx Only	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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