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**ENGLISH**

**Indication:**

Polyester abrasive strips are instruments used in dental procedures for finishing, polishing, and interproximal wear. They allow for controlled removal of excess restorative material and adjustment of the contact point between adjacent teeth.

**Intended user:**

Duly qualified dental professionals.

**Target population:**

The product is intended for use in dental patients undergoing restorative or orthodontic procedures (interproximal wear) as indicated by the dentist.

**Fundamentals of its operation and action:**

Polyester abrasive strips act by abrasive mechanical action, promoting the controlled removal of excess restorative material or tooth enamel through abrasion between the active surface of the strip and the area to be worn. This action allows for interproximal wear, finishing, and polishing of restorations or dental surfaces.

**Contraindications:**

- Do not use in patients with known hypersensitivity to metals, especially stainless steel.
- Do not use in areas with severe inflammation or active bleeding.
- Do not use the product that is damaged, deformed or has lost its abrasive coating.

**Instructions for use:**

Polyester abrasive strips feature a neutral center (non-abrasive) that is utilized to insert the sandpaper into the proximal faces without wearing down the contact area.

- Hold the sandpaper with the pulps of your thumb and index fingers.
- Insert the sandpaper with the neutral center into the interproximals of the teeth.
- Make movements from the buccal to palate/lingual and vice versa, from the cervical region towards the occlusal surface and vice versa.
- Finishing and polishing should start on the darker side of the strip (medium abrasiveness) and finish on the lighter side (fine abrasiveness).

**Form of presentation:**

- Non-sterile product.
- One-time use.
- Pack containing 50, 100 or 150 units.
- 3 models: 2.5 mm, 4.0 mm or 6.0 mm.

**Cleaning Procedure:**

Not applicable. The product is for single use, so it should not be cleaned or reused.

**Sterilization Procedure:**

Not applicable. The product is supplied non-sterile and should not be sterilized. Dispose of dental waste after use in accordance with local regulations.

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REV.01

**Warnings and Precautions:**

Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.

- Goggles should be worn to protect against ejected particles.
- Wear PPE, including gloves, masks, safety goggles, and long-sleeved lab coats when handling the Strip, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

**Conservation and storage:**

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The operator must be informed of the content. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the instrument from shipping to delivery to the customer. Packets should not be left to the action of time.

**Discard:**

- The disposal of the material must follow the relevant local legislation.

**Shelf Life:**

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

**Service life:**

Not applicable. The product is for single use only and should be discarded after use.

**Additional information:**

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (Qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: [www.microdont.com.br](http://www.microdont.com.br), the most up-to-date version of the instructions for use corresponding to the product.

**Composition:**

Polyester and Aluminum Oxide.

**Symbology:**

| Symbol                                                                              | Symbol Title   | Description                                                                                                     |
|-------------------------------------------------------------------------------------|----------------|-----------------------------------------------------------------------------------------------------------------|
|  | Use-by date    | Indicates the date after which the medical device should not be used.                                           |
|  | Non-sterile    | Indicates a medical device that has not undergone a sterilization process.                                      |
|  | Do not reuse   | Indicates a medical device intended for a single use, or for use on a single patient during a single procedure. |
|  | Medical device | Indicates that it is a medical device.                                                                          |

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REV.01

|                                                                                   |                                                     |                                                                                                                  |
|-----------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
|  | CE Mark                                             | It indicates that the medical device is in compliance with Regulation (EU) 2017/745 (MDR).                       |
|  | Consult instructions for use                        | Indicates the user's need to refer to the instructions for use.                                                  |
|  | Date of manufacture                                 | Indicates the date the medical device was manufactured.                                                          |
|  | Manufacturer                                        | Indicates the manufacturer of the medical device as defined in Regulation (EU) 2017/745 (MDR).                   |
|  | Lot Code                                            | Indicates the manufacturer's batch code so that the batch or lot can be identified.                              |
|  | Catalog Number                                      | Indicates the manufacturer's reference number so that the medical device can be identified.                      |
|  | Authorised representative in the European Community | Indicates the authorized representative in the European Community.                                               |
|  | For professional use only                           | It indicates that the medical device must be used exclusively by a qualified health professional.                |
|  | Prescription only                                   | Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation. |
|  | Unique Device Identification (UDI)                  | Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).                              |
|  | QR Code – Access to the electronic IFU              | QR Code to access the electronic version of the product's Instructions for Use (eIFU).                           |

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