CUNHA PLÁSTICA
PLASTIC DENTAL WEDGE
CUÑA DENTAL PLÁSTICA

REV.01

ENGLISH

Indication:

Anatomical wedges assist in restorative procedures of class II cavities, that is, those that involve the proximal surface of the teeth, have the mechanical function of adapting the matrix and/or moving adjacent teeth apart. Because they are anatomical, they adapt to the papilla, greatly minimizing the chance of damage to the gum tissue.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures involving proximal surfaces of the teeth, in which it is necessary to stabilize and properly adapt the dental matrices.

Fundamentals of its operation and action:

The acute shape of the anatomical wedge facilitates its introduction into the embrasure, fixing the matrix and helping to avoid excess restorative material.

Anatomical Wedge Sizes:

White: Large

Light Blue: Medium

Dark Blue: Small

Contraindication:

When there is a proven allergy to any of its components, the product should not be used.

Instructions for Use:

- Select the size wedge compatible with the clinical situation.
- Insert the wedge into the interproximal space, usually from lingual to buccal, due to the larger size of the embrasure.
- The anatomical wedges resemble an isosceles triangle, so when inserted, their faces of the same size should face the proximal surfaces of the teeth and the face of the wedge of a different size facing the gum, this face has anatomy to be accommodated over the papilla, without causing damage to it.
- Check if the matrix has been correctly adapted, or if there will be a need to replace it with a larger wedge.
- Perform the restoration using the desired matrix system.
- Remove the wedge after completion of the restoration procedure.

CUNHA PLÁSTICA
PLASTIC DENTAL WEDGE
CUÑA DENTAL PLÁSTICA

REV.01

Cleaning Procedure:

Not applicable. The product is for single use only and should not be cleaned or reprocessed after use.

Sterilization Requirements

As the product is supplied non-sterile and intended for single use, it must not be sterilized, autoclaved or reprocessed, and must be discarded after completion of the procedure. Reuse is prohibited.

Warnings and Precautions:

Only properly trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.

- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- It is recommended that the use of the product is always accompanied by absolute isolation of the operative field, since small objects can be swallowed or aspirated accidentally.
- Non-sterile product – do not sterilize or reprocess.
- Do not use the product if the packaging is damaged.
- Single use – discard after use.
- Its reuse may cause cross-contamination, forbidden reuse.
- Follow current legislation for the disposal of dental products.
- The product has a pointed end and should be handled in such a way as not to cause injury to the patient.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

CUNHA PLÁSTICA
PLASTIC DENTAL WEDGE
CUÑA DENTAL PLÁSTICA

REV.01

Composition:

Material: Plastic

Symbology:

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	CE Mark	It indicates that the medical device complies with Regulation (EU) 2017/745 (MDR).
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product as defined in Directives 90/385/EEC, 93/42/EEC and 98/79/EC of the European Community.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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