CUNHA PROTETORA
PROTECTIVE DENTAL WEDGE
CUÑA DENTAL PROTECTORA

REV.01

ENGLISH

Indication:

Protective wedges, as their name suggests, are indicated to protect adjacent teeth in restorations, such as in class II clinical situations, it reduces the risk of damage to the neighboring tooth and at the same time performs the function of an anatomical wedge.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the use of protective wedges to protect adjacent dental structures, as indicated by the dentist.

Fundamentals of its operation and action:

The sharp shape of the protective wedge facilitates its introduction into the embrasures, and the metal part that runs along the entire body of the wedge protects the adjacent tooth.

Wedge sizes:

Lilac: Large

Orange: Medium

Green: Small

Contraindication:

When there is a proven allergy to any of its components, the product should not be used.

Instructions for Use:

- Select the size wedge compatible with the clinical situation.
- Insert the wedge into the interproximal space, usually from lingual to vestibular, due to the larger size of the niche.
- The protective wedges are shaped like an isosceles triangle. That way, when inserted, its equal-sized faces should face the proximal surfaces of the teeth, while the different-sized face should face the gingiva.
- Check if the metal part of the wedge has been correctly adapted, or if there will be a need to replace it with a larger wedge.
- Perform cavity preparation with the protective wedge in position, to prevent the diamond tip or drill from touching the healthy tooth surface of the adjacent tooth.
- After preparation, remove the protective wedge and perform the restoration using the desired matrix system.

Cleaning Procedure:

CUNHA PROTETORA
PROTECTIVE DENTAL WEDGE
CUÑA DENTAL PROTECTORA

REV.01

Not applicable. The product is for single use only and should not be cleaned or reprocessed after use.

Sterilization Requirements

As it is a disposable device, after use, it must be discarded, and its reuse is unfeasible.

Warnings and Precautions:

- Only properly trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.
- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- It is recommended that the use of the product is always accompanied by absolute isolation of the operative field, since loose objects can be swallowed or aspirated accidentally.
- Non-sterile product, autoclaved prior to use.
- The product has a pointed end and should be handled in such a way as not to cause injury to the patient. Professional use.
- Do not use the product if the packaging is damaged.
- Single use – discard after use.
- Its reuse may cause cross-contamination, forbidden to reprocess.
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified on the product label, to be informed of any occurrences with the manufacturer.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and should be discarded after use.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Stainless steel and plastic

CUNHA PROTETORA
PROTECTIVE DENTAL WEDGE
CUÑA DENTAL PROTECTORA

REV.01

Symbology:

Symbol	Symbol Title	Description
	Expiry date	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates a health product intended for a single use, or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	CE Mark	It indicates that the medical device complies with Regulation (EU) 2017/745 (MDR).
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product as defined in Regulation (EU) 2017/745 (MDR).
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identifier (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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