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ENGLISH

Indication:

The clamp (also known as the die fixation ring) is an indispensable component for the use of sectional dies, being fundamental for the success of dental restorations. When used correctly in class II restorations, which involve the proximal surfaces of the dental element, such as mesio-occlusal-distal (MOD), mesio-occlusal (OM) and occlusal-distal (DO), the clamp has the main function of fixing the sectional matrix to the proximal wall of the tooth to be restored. In this way, the clamp ensures the stabilization and proper conformation of the matrix, allowing the precise performance of the restorative procedure.

3R Indústria e Comércio Ltda is not responsible when the product is used outside the application for which it was designed.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures that involve the interproximal aspect of the teeth, such as Class II restorations, in which it is necessary to use sectional matrices and staples for their fixation, according to the clinical indication of the dental surgeon.

Fundamentals of its operation and action:

Mechanical: The clamp, along with sectional dies and wedges play a crucial role in dental restoration, aiding in the conformation and stabilization of the die.

Clip types:

Clamps that do not require forceps (Model A): Available in a universal size that meets clinical needs; it has a stem for finger grasping.

Clamps that require forceps (Model B): Available in a universal size that meets clinical needs; it requires clamp-holding forceps (the same type used to grip clamps for rubber dam isolation).

Model A

Model B

Contraindications:

- It should not be used in patients with allergy or hypersensitivity to any component material of the product.
- It is not suitable for teeth with fragile, compromised structures or with extensive loss of dental tissue, where fixation can cause additional damage.
- Avoid use in patients with marked tooth mobility or in deciduous teeth close to exfoliation, due to the risk of displacement or discomfort.

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REV.01

- Do not use in clinical situations where there is no adequate access for safe positioning of the staple, which can cause gingival lesions or damage to the restoration.

Instructions for use:

In the case of a restoration that involves both proximal faces (e.g., MOD), the clamp secures the sectional matrix and makes it possible to perform the procedure on one face at a time. This is very different from the steel die band with the Tofflemire die holder, as it improves the field of view of cavity preparation and makes it easier to perform the job with greater precision.

- Select and insert the protective wedge, to wear the face of the tooth in question more safely;
- Once the preparation is finished, remove the protective wedge with the help of tweezers;
- If necessary, perform pre-wedge or slight wear with a sandpaper strip on the prepared tooth to allow accommodation of the die without denting it;
- Select the sectional matrix that best suits the tooth to be restored;
- If pre-minting has been performed, remove the wedge to allow the die to fit and then reposition it according to the space in question;
- Choose the clamp that is compatible with the clinical situation at hand.
- With the clamp, fix the matrix on the tooth surface. In extreme cases, the conformation of the die can be adjusted with the help of a hone or explorer probe, but be careful at this stage not to dent the die that already has an ideal shape;
- Once the procedure is completed, remove the clamp, then the wedge and finally the die.

Manual Cleaning Procedure

If manual cleaning is the only option available, proceed as follows:

- Use an area specifically reserved for cleaning.
- Rinse the device under running water and keep it immersed in a bath of neutral pH cleaning solution (enzymatic detergent), following the instructions of the agent manufacturer.
- After soaking and keeping it immersed, brush well away from the body using the neutral cleaning agent.
- Care should be taken to avoid spreading contaminants by spraying during the brushing process.
- Do not use wire brushes.
- Special care should be taken to clean crevices and other hard-to-reach areas. Visually inspect to confirm debris removal.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.
- Repeat the cycle if necessary.

Ultrasonic Cleaning Procedure

For cleaning by means of ultrasound, it is recommended to proceed as follows:

- Prepare a pH-neutral cleaning solution (enzymatic detergent), following the agent manufacturer's instructions for the correct concentration, exposure time, temperature, and water quality.
- Place in the proper location of the ultrasonic equipment, following the equipment manufacturer's instructions for use.
- Completely submerge the device in the cleaning solution and program the equipment for at least fifteen (15) minutes.
- Rinse thoroughly under running water until visibly clean.
- Dry using a clean compressed air, or paper towel.

Sterilization Requirements

- Sterilization by pressure saturated steam autoclave at 121°C for 20 minutes is recommended, or sterilization parameters should be followed in accordance with the relevant policies/regulations of the country where the device is being used. It is the

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REV.01

professional's full responsibility to keep the autoclave equipment in ideal conditions of use to ensure the correct sterilization of the product.

- For sterilization, remove the product from its original packaging, pack in material suitable for sterilization, and sterilize in autoclave (humid heat, high temperature and pressure) according to the equipment's instructions.

Warnings and Precautions:

- Only dental surgeons should use this product, as they have the scientific information necessary for its proper use.
- Before use, inspect whether the clamp is in perfect condition for use;
- Non-sterile product.
- Sterilize the product before use.
- During the operation, the professional must use all personal protective equipment, including protective goggles.
- Wear protective goggles on the patient.
- LOT Number: It is important to keep the lot number, identified on the product label, to be informed in any occurrences with the manufacturer.

Discard:

- The disposal of the material must follow the relevant local legislation.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Life:

The shelf life of the product may vary according to the conditions of use, frequency of sterilization, storage and maintenance. It is recommended that the product be regularly inspected for signs of wear, deformation or corrosion, and replaced whenever it presents any alteration that may compromise its performance or safety.

Conservation and storage:

- Keep the packaging in a dry place, away from light and away from the action of bad weather.
- Keep the products in their original packaging until the moment of use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you no later than 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

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REV.01

	Expiry date	Indicates the date after which the medical device is not to be used.
	Non-sterile	Indicates a medical device that has not undergone a sterilization process.
	Keep dry	Indicates a medical device that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the medical device was manufactured.
	Keep out of sunlight	Indicates a medical device that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the medical device.
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the medical device can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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