ENGLISH**Indication:**

The disposable brush applicators are indicated for the controlled, precise and uniform application of dental substances and materials during clinical procedures. They may be used for the application and spreading of caries and dental plaque disclosing solutions, adhesive systems, sealants, haemostatic agents and fluorides, among other materials compatible with the clinical technique. The brush tip allows the deposition and distribution of the materials on the tooth surfaces and other areas of interest, contributing to greater control and efficiency during application. The flexible, bendable shaft provides better access to different regions of the oral cavity, favouring precision and ease of handling during the procedure.

Intended user:

Duly qualified dental professionals.

Target population:

Patients undergoing dental procedures involving the topical application of solutions, agents or fluid materials, regardless of age or clinical condition, according to the indication and judgement of the responsible professional.

Mechanism of action:

The disposable brush applicator acts by retaining and transferring dental substances and materials through its nylon bristles, allowing their controlled deposition and distribution on the tooth surface or other areas of interest during the clinical procedure. The flexibility and arrangement of the bristles enable uniform application of the material, while the shaft provides support and control during handling of the applicator.

Contraindications:

- Do not use in patients with known hypersensitivity to the components of the device.
- Do not use with highly caustic substances or strong solvents that may compromise the integrity of the fibres or of the shaft.
- Do not use for the application of highly viscous or very thick materials that may deform or damage the active part of the brush applicators.

Instructions for use:

- Remove the applicator from the packaging using an appropriate technique to avoid contamination of the bristles.
- If necessary, adjust the bendable shaft of the applicator according to the region of the oral cavity to be accessed.
- Insert the active part into the solution to be applied.
- Position the bristles on the surface or region to be treated and make gentle, controlled movements to deposit and evenly distribute the material.

Presentation:

- Non-sterile product.
- Single use.
- Package containing 100 units.
- Available in regular size only.

Cleaning procedure:

Not applicable. The product is for single use and must therefore not be cleaned or reused.

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REV.00

Sterilization procedure:

Not applicable. The product is supplied non-sterile and must not be sterilized. Discard after use in accordance with local dental waste regulations.

Warnings and Precautions:

- Product for professional use only.
- Keep out of the reach of children.
- Do not use the product if the packaging is breached, damaged or previously opened.
- Do not expose the device to high temperatures.
- During the procedure the professional must use all personal protective equipment, including protective eyewear.
- Ensure that the patient also wears protective eyewear.
- LOT number: it is important to keep the lot number, shown on the product label, so that it can be reported to the manufacturer in the event of any occurrence.

Conservation and storage:

- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The product must be stored in a well-ventilated, dry place protected from light, so as to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.

Disposal:

- After use, discard the applicator in an appropriate container, in accordance with the classification adopted in the Health Service Waste Management Plan (PGRSS) of the generating facility and with ANVISA Resolution RDC No. 222/2018, or in accordance with the local legislation applicable in the country of use.

Shelf life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the indicated conditions, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use and must be discarded after use.

Other additional information:

- If necessary, 3R provides the instructions for use in printed format, at no additional cost, upon request by e-mail (qualidade@3R.ind.br) or by telephone (+55 11-5525-0590 / +55 0800-773-5525), and we will reply within a maximum of 7 business days from the request.
- Always check the website www.microdont.com.br for the most up-to-date version of the instructions for use corresponding to the product.
- In the event of adverse events or any incidents related to the use of the product, immediately notify the manufacturer by e-mail: qualidade@3R.ind.br.

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REV.00

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	CE Mark	Indicates that the medical device complies with Regulation (EU) 2017/745 (MDR).
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the medical device manufacturer, as defined in Regulation (EU) 2017/745 (MDR).
	Lot Code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	Authorised representative in the European Community	Indicates the authorized representative in the European Community
	Prescription only	Sale and use restricted to qualified professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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