ENGLISH

Indication:

Wooden wedges are indicated in restorative procedures with the mechanical function of adapting the matrix and/or moving adjacent teeth apart.

Intended user:

Duly qualified dental professionals.

Target population:

The product is intended for use in dental patients undergoing restorative procedures involving proximal surfaces of the teeth, in which it is necessary to stabilize and properly adapt the dental matrices.

Fundamentals of its operation and action:

The sharp shape of the wedge facilitates its introduction into the embrasure, fixing the matrix and helping to avoid excess restorative material.

Contraindication:

- When there is a proven allergy to any of its components, the product should not be used.
- It is recommended that the use of the product is always accompanied by absolute isolation of the operative field, since loose objects can be swallowed or aspirated accidentally.
- The product has a pointed end and should be handled in such a way as not to cause injury to the patient.
- Professional use.

Instructions for Use:

- Select the desired wedge.
- Insert the wedge into the interproximal space, usually from lingual to vestibular, due to the larger size of the embrasure.
- Wooden wedges resemble an isosceles triangle, so when inserted their faces of the same size should face the proximal surfaces of the teeth and the face of the wedge of a different size facing the gums.
- Check whether the matrix has been correctly adapted, or whether it will be necessary to replace it with a larger wedge.
- Perform the restoration using the desired matrix system.
- Remove the wedge after completion of the restoration procedure.

Form of presentation

- Non-sterile product;
- One-time use

N°1 – 12mm X 2.2mm – green

N°2 – 13mm X 2.2mm – light blue

N°3 – 14mm X 2.5mm – lilac

N°4 – 15mm X 2.5mm – orange

N°5 – 16mm X 2.5mm – natural

Warnings and Precautions:

CUNHAS ODONTOLÓGICAS DE MADEIRA
WOODEN DENTAL WEDGES
CUÑAS DENTALES DE MADERA

REV.01

- Only properly trained and qualified dental professionals should use this product, as they have the scientific information necessary for its proper use.
- Wear PPE, including gloves, masks, safety glasses, and long-sleeved lab coats when handling the product, as is done with all other dental instruments. Wear protective glasses on the patient as well.
- Do not use the product if the packaging is damaged, opened or wet.
- Manufacturer recommends single use – discard after use.
- Its reuse can cause cross-contamination.
- Follow current legislation for the disposal of dental products.
- BATCH NUMBER: It is important to keep the batch number, identified in the product title, to be informed of any occurrences with the manufacturer.

Conservation and storage:

- The product must be stored in an airy, dry place protected from light, in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity.
- Keep the product in its original packaging until the moment of use. Inspect the packaging before use.
- The carrier must be informed about the contents. The precautions to be taken in transport are demonstrated on the packaging, thus ensuring the protection of the product from shipping to delivery to the customer. Packets should not be left to the action of time.

Discard:

The disposal of the material must follow the relevant local legislation.

Shelf Life:

Not applicable. The product is manufactured from non-perishable materials and supplied non-sterile. When stored under the conditions indicated, it retains its properties indefinitely.

Service life:

Not applicable. The product is for single use only and must be discarded after use.

Additional information:

- If necessary, 3R provides instructions for use in printed form, at no additional charge, upon request by email (qualidade@3r.ind.br) or by phone: (+55 11-5525-0590 / +55 0800-773-5525), and we will get back to you within a maximum of 7 business days from the request.
- Always check the website: www.microdont.com.br, the most up-to-date version of the instructions for use corresponding to the product.

Composition:

Colored dental wedges:

Wood (structural function)

Non-toxic reactive dyes (pigmentation function) 1%

Symbology:

CUNHAS ODONTOLÓGICAS DE MADEIRA
WOODEN DENTAL WEDGES
CUÑAS DENTALES DE MADERA

REV.01

Symbol	Symbol Title	Description
	Validity	Indicates the date after which the health product should not be used.
	Non-sterile	Indicates a health product that has not undergone a sterilization process.
	Do not reuse	Indicates an intended health product for use or for use on a single patient during a single procedure.
	Keep dry	Indicates a health product that needs to be protected against moisture.
	Medical device	Indicates that it is a medical device.
	Recyclable packaging	The packaging is recyclable.
	See instructions for use	Indicates the user's need to refer to the instructions for use.
	Date of manufacture	Indicates the date the health product was manufactured.
	Keep out of sunlight	Indicates a health product that needs protection from light sources.
	Manufacturer	Indicates the manufacturer of the health product.
	Lot Code	Indicates the manufacturer's batch code so that the batch or batch can be identified.
	Catalog Number	Indicates the manufacturer's reference number so that the health product can be identified.
	For professional use only	It indicates that the medical device must be used exclusively by a qualified health professional.
Rx Only	Prescription only	Sale and use restricted to licensed professionals or by prescription, in accordance with applicable legislation.
	Unique Device Identification (UDI)	Indicates that the adjacent code corresponds to the Unique Device Identifier (UDI).
	QR Code – Access to the electronic IFU	QR Code to access the electronic version of the product's Instructions for Use (eIFU).

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